## grant



## Refrigerated and Heated Circulators and Chillers

LT ecocool Energy efficient refrigerated heated circulating baths

Optima range Refrigerated heating circulating bath combinations

RC series Recirculating chillers

# Refrigerated and heating circulators and chillers

Cost-effective and efficient multi-purpose systems for cooling applications. The circulating product portfolio offers a diverse choice to meet a range of needs and budgets. We offer entry-level products for standard use to advanced-level products for more demanding requirements or opt for a product solution that is customised to your individual needs.

Energy efficient models - delivers powerful cooling and significant running cost savings

Powerful precision cooling - for use in open-loop or closed-loop format

User-friendly and intuitive design - handy features, effortless maintenance and compact design

Robust, durable construction - for longevity, reliability and long-term low cost of ownership

A comprehensive range - multiple combinations or custom solutions to meet your specific needs

Added protection - industry leading warranty options up to four years





#### LT ecocool range

Grant R4

### Operating Temperature

The LT ecocool refrigeration range offers accurate temperature control from -25°C to 150°C and is available in two models.

The R4 refrigeration unit can be combined with any of Grant's Optima<sup>™</sup> heated circulators to offer a temperature range of up to -30°C to 100°C as standard. Custom units can be designed for wider temperature ranges.

## Six points to consider when choosing your system

#### Do you need to immerse samples within a tank?

Consider the working area required. The table on page 27 shows the dimensions of the top opening and the min/max fluid depths.

#### Cooling power required at a given temperature

For example, if your operating temperature is 0°C and you need 500W cooling power, you need the R4 refrigeration unit with any of the controllers. Alternatively to calculate the power required use the following formula:

 $W = V \times T \times K / 60 \times t$  (mins)

#### Cool-down time required to reach that temperature

Calculate the cool-down time required according to the following formula and refer to the cool down curved for individual performance.

 $W (mins) = V \times T \times K/60 \times t (mins)$ 

W = average cooling power	Water	K = 4180
V = total system fluid volume L	50/50 water/glycol	K = 3800
T = temperature difference °C	Alcohol	K = 2100
K = fluid heat capacity (j/L°C)	Silicone oil	K = 1800

## Do you need to control the temperature of or remove the heat from an external device?

1. Consider the pump requirement. Fluid flow rate is critical in order to maintain adequate exchange of heat within the external system. Flow rate is dependent on the restrictions within the system. Factors which cause a pressure drop are height, length, pipe bore and the number and angle of bends within the system. To maintain sufficient flow in a highly restricted system, a high pressure pump is required. The integral pumps in the Optima<sup>™</sup> and LT ecocool series thermostats are satisfactory for most laboratory applications; for more powerful pump requirements select either of the Grant accessory vertical turbine pumps (VTP).

2. Consider whether you need to control the temperature within the external apparatus. For external temperature control choose the TX150, TXF200 or LT ecocool 150 controllers and an external temperature probe.

#### Do you require temperature ramping?

If yes, choose the TX150, TXF200 or LT ecocool controller and Labwise® accessory software. For refrigeration on/off choose refrigeration units LT ecocool 150 or R4.

#### What other features do you require?

Compare the wide range of features offered by the four  $Optima^{TM}$  series or LT ecocool 150 controllers and select the controller that meets your needs.

#### Need more help?

If you need help choosing the correct system, please contact: salesdesk@grantinstruments.com or call +44 (0) 1763 260 811.



## LT ecocool Energy saving refrigerated and heated circulating baths

The range of innovative, eco-friendly, refrigerated heated circulating baths reduce operating costings and help to protect the environment by achieving energy savings of up to 80%\*. A choice of two models, both supplied assembled as ready-to-use kits, complete with accessory tubing\*, clips and connectors as standard.

Choice of two models, temperature range -25 to 150°C\* (vary on model choice)

#### Industry leading 4 year warranty with online registration

Active cooling throughout the full temperature range

Energy savings of up to 80% \*compared to standard compressor units



\* Temperature range of tubing supplied: -40°C to 100°C (can be length as required). Supplied tubing 2 x 1.5m ID 9mm ø

### Applications

- Pharmaceutical mini pilot plant reactors
- Education immersing small samples, photometry, chromatography systems
- Industrial QC testing, sample preparation, general cooling, reaction chemistry, temperature control, semi-conductor manufacturing, rheometry
- Food refractometry
- Life-science electrophoresis cooling

## LT ecocool refrigeration range Technical specifications

		LT ecocool 100	LT ecocool 150	
Dimensions	hxdxw mm	640 x 4	30 x 245	
Capacity	L	7	7	
Temperature range	°C	-20 to 100	-25 to 150	
Stability	±°C	0.05	0.02	
Flow rate (max)	L/min	17	14-22 (adjustable)	
Pump pressure (max)	mbar	250	530	
Working area	d x w mm	118 >	< 154	
Min/max fluid level	mm	85/145	130/145	
Calibration points		2	5	
Cooling power (mean)	@20°C W	227	415	
	@0°C W	190	227	
	@-10°C W	115	117	
	@-20°C W	41	71	
Programmes		-	1 x 30 segments Labwise® required	
Communication interface		-	USB	
Temperature probe socket		-	6 pin mini DIN	
Display		4 digit LED	Full colour QVGA TFT	
Languages		-	EN, FR, DE, IT, ES	
Timer		1 minute to 99 h	ours 59 minutes	
Temperature presets			3	
Alarms		High	High and low	
Electrical supply (max) A	220-240V	12 (50 or 60 Hz versions available)	12 (50/60Hz)	
	110-120V	18.5 (50/60Hz)	18.5 (50/60Hz)	
Safety		Adjustable over te	mperature cut-out	
Ready to use kits		Assembled and supplied with standard tubing, insulation, clips and connectors		
Weight	kg	2	9	

## Fluids

We recommend the following fluids for use in Grant baths:

-50°C to 50°C:	Silicone oil - Iow viscosity Bayer silicone M3			
-30°C to 70°C:	50% water, 50% antifreeze - inhibited ethylene glycol			
0°C to 30°C:	80% water, 20% antifreeze - inhibited ethylene glycol			
5°C to 99.9°C:	Water - do not use to boil water			
70°C to 150°C:	Silicone fluid (viscosity ~20cs, flash point ≥230°C, fire point ≥280°C)			
Always read the manual and warnings when choosing a fluid.				

## Optima<sup>™</sup> Refrigerated baths and circulator range

High-performance refrigeration unit easily combined with any of our four Optima<sup>™</sup> heated circulators. Offers flexibility and delivers outstanding temperature performance for routine and more sophisticated applications requiring accurate temperature control in the range of -30°C to 100°C. Also available as a kit, Grant offer the LTC4 (TX150-R4) with the heated circulator, refrigeration unit and insulated tubing\* offering a complete ready-to-use system.

Choice of two base refrigeration units and four heated circulators, temperature range -30°C to 100°C\*\* (vary on model choice)

Stability: Up to ±0.01°C

No spill drain valve located on the front of the unit

Safe - water freeze protection thermostat and 27 bar high pressure switch

Three pre-set programs

3 years warranty, 4 years with the LTC4

#### Custom units are available



\* Temperature range of tubing supplied: -40°C to 100°C (can be cut to length as required).

\*\* Lower temperatures are available. Contact +44 (0) 1763 260 811 or email salesdesk@grantinstruments.com to find out more.

## Grant R series base refrigeration units

Technical specifications

		R4	LTC4
Dimensions	hxdxw mm	550 x 515 x 393	755 x 515 x 393
Capacity	L	20	20
Temperature range (T100)	°C	0-100	-
Temperature range (TC120)	°C	-20 to 100	-
Temperature range (TX150)	°C	-30 to 100	-30 to 100
Temperature range (TXF200)	°C	-30 to 100	-
Refrigerant		R134a	R134a
Working area	d x w mm	230 x 305	230 x 305
Min/Max fluid level	mm	85/140	85/140
Cooling power (typical)	@20°C W	900	900
	@0°C W	500	500
	@-10°C W	300	300
	@-20°C W	180	180
	@-30°C W	40	40
	@-40°C W	-	-
	@-47°C W	-	-
Electrical power (max) W	120V	780 (50-60Hz)	2280 (50-60Hz)
	230V	850 (50Hz)	2850 (50Hz)
Relay control*		•	•
Weight	kg	40.6	42.9

\* relay to enable switching off the refrigeration system in a program

## Applications

• University research/teaching - temperature control of external equipment including: spectrophotometers & refractometers. Circulation of temperature control fluid to jacketed vessels, cooling crystallisation vessels

• Industrial laboratories - temperature probe calibration, product testing, product QC, temperature control of external equipment.

## Grant Optima<sup>TM</sup> heated circulators Technical specification

	Grant Op	tima Heated circulate	ors and Immersion t	hermostat		
	T100	TC120	TX150	TXF200		
Dimensions h x d x w	333 x 172 x 120	333 x 172 x 141	342 x 1	72 x 141		
Stability (DIN 12876) water @10°C	⊧°C	C	).]			
Setting resolution	°C	D.1	0.1 (0.01 wit	h Labwise®)		
Programs		-	1 x 30 segments Labwise® required	10 x 100 segments		
Safety over tempera	ure Fixed	Fixed Adjustable cut-out				
Alarms (can be configured to switch a relay)	-	- High (no relay)		High and low		
Language capability		-		EN, FR, DE, IT, ES		
Height above tank rim	nm	200				
Depth below tank rim	nm	13	35			
Display	4 dig	it LED	Full colour	QVGA TFT		
Timer	-	lm	inute to 99 hours, 59 minu	utes		
Calibration points		2		5		
Communication interface		-		USB, RS232, remote temperature probe		
Heater power W 120V/2	30V 144C	1440/1290		/1840		
Electrical power (50/60Hz) W 120V/2	30V 150C	1500/1400 1500/2000		/2000		
Weight	kg 2.1	2.3	2.6	2.6		

For more information on the Grant Optima heated circulators, please see page 1.8.

#### Options and accessories

Labwise® PC	software		
Allows two-way communication for status display, programming and data capture (see page 3.1 for more information) USB/RS232 cables provided.			Compatible with TX150, TXF200 and LT ecocool 150 models.
External prob	oes for monitori	ng and control	lling temperature of remote loads
TXPEP flexible pl	astic probe, 3m cab	le (Din plug)	Compatible with TX150 and TXF200 models. Compatible with LT ecocool 150 models
TXSEP stainless s	teel probe, 3m cabl	e (Din Plug)	Compatible with TX150 and TXF200 models. Compatible with LT ecocool 150 models
Vertical turbi	ne pumps*		
	bact design. Supplie d special lid for fittin		Required only where application demands a higher pressure than that delivered by the internal pump to maintain flow.
VTP 1 Max pressure Max. flow	1000 mbar 9 L/min	230V 50Hz	Note: The optional VTP pumps will transfer additional heat to the baths and reduce the net cooling power of the refrigeration unit. The above figures must be taken into consideration when choosing the refrigeration unit. When order a VTP pump, please specify which refrigeration base unit it is to be used with.
VTP 2 Max pressure Max. flow	1650 mbar 12 L/min	230V 50Hz	Note: Other sizes of heat exchange coil can be made to your specification, contact Grant for further information.
Heat exchan	ge coil		
CW5 Other sizes of heat exchange coil can be made to your specification, contact us for further information.			Temperature range: 2°C above the temperature of the coolant Coil Øxl: 77 x 55mm Pipe bore inlet/outlet: 7mm
Hose kits			
HOSE100 General purpose tubing and insulation kit: -40°C to 100°C			Tube kit 2 x 2m, assembled with Optima pump outlet plate and simple tube clips,
HOSE200 High to kit: -50°C to 200	emperature tubing ; P°C	and insulation	no tools required.

## High pressure pumps Optional

		VTP p	Heat Exchange Coil	
		VTP-1LT	VTP2-LT	CW5
Maximum pressure	water mbar	1000	1650	-
Maximum flow	water L/min	9	12	-
Pipe bore	inlet/outlet mm	12	2.7	7
Electrical connection		10 am	-	
Power consumption	W	30	40	-
Power output to fluid	W	15	22	-
Safety		Therm	al fuse	-
Temperature range	°C		2°C above coolant temperature	
Coil Øxl	mm		_	77 x 55

## Grant refrigeration units Optional

				1925					
			LT ecocool 100	LT ecocool 150	T100-R4	TC120-R4	TX150-R4	TXF-200-R4	LTC4
Labwise® Software (see section 3 for furthe	r informatic	on)	-	•	-	-	•	•	•
CW5 Heat exchange co	oil		•	•	•	•	•	•	٠
IQOQ Documentation			IQOQ LT ecocool 100	IQOQ LT ecocool 150	IQOQ T100 + IQOQ R4	IQOQ TC120 + IQOQ R4	IQOQ TX150 + IQOQ R4	IQOQ TXF200 + IQOQ R4	IQOQ LTC4
PQ Documentation			PQ LT ecocool 100	PQ LT ecocool 150	PQ T100 + PQ R4	PQ TC120 + PQ R4	PQ TX150 + PQ R4	PQ TXF200 + PQ R4	PQ LTC4
Extended warranty 1	year	EWC1	•	٠	•	•	٠	•	٠
Extended warranty 2	years	EWC2	-	-	•	•	•	•	-
Temperature pro	bes - 3n	n cable							
TXPEP Plastic probe					-		•	•	•
TXSEP Stainless steel	probe				-		•	•	•
PEP Plastic probe			-	•			_		
SEP Stainless steel pr	robe		_ · ·						
Pumps - optional									
VTP1-LT Maximum pressure Maximum flow	1000 mb 9 L/min	ar					_		
VTP2-LT Maximum pressure Maximum flow	1650 mba 12 L/min		•		-				
VTP1-PLR4 Maximum pressure Maximum flow	1000 mk 9 L/min	bar	-	-			•		
VTP2-PLR4 Maximum pressure Maximum flow	1650 mba 12 L/min		-	-			•		

## Pumps

P-M6	ofto	Replacement plastic pump inlet/outlet connector. Fits tubing 9mm inner dia. Temperature range -50°C to 200°C.
P-M11	0000	Replacement plastic pump inlet/outlet connector. Fits tubing 15mm inner dia. Temperature range -50°C to 200°C.
M-SR4		Metal pump inlet/outlet connector, dual seal super rapid 4mm. Fits semi rigid tubing 4mm outer dia. Temperature range -20°C to 100°C.
M-SR6		Metal pump inlet/outlet connector, dual seal super rapid 6mm. Fits semi rigid tubing 4mm outer dia. Temperature range -20°C to 100°C.
M-SR8		Metal pump inlet/outlet connector, dual seal super rapid 8mm. Fits semi rigid tubing 4mm outer dia. Temperature range -20°C to 100°C.
M-HB7	all the	Metal pump inlet/outlet connector, tube barb 7mm. Fits flexible tubing 7mm inner dia. Temperature range -50°C to 200°C.
M-HB9		Metal pump inlet/outlet connector, tube barb 9mm. Fits flexible tubing 9mm inner dia. Temperature range -40°C to 120°C.
M-HB12	ale les	Metal pump inlet/outlet connector, tube barb 12mm. Fits flexible tubing 12mm inner dia. Temperature range -40°C to 120°C.
M-UC	40000	Metal pump inlet/outlet plate, 1/4" BSP/G1/4 female. Temperature range -50°C to 200°C.
HOSE100		General purpose tubing and insulation kit, includes 2 x 2m general purpose insulated tubing -40 to 100°C, assembled with LT ecocool/Optima <sup>™</sup> pump outlet plate and simple tube clips, no tools required. Can be cut to length. 10mm ID, 14mm OD.
HOSE200		General purpose tubing and insulation kit, includes 2 x 2m general purpose insulated tubing -50 to 200°C, assembled with LT ecocool/Optima™ pump outlet plate and simple tube clips, no tools required. Can be cut to length. 8mm ID, 11mm OD.

## RC Series Recirculating chillers

A choice of two robust recirculating chillers delivering a constant flow of temperature-controlled fluid to provide powerful, regulated cooling at -10°C for many types of industrial machinery and scientific apparatus. Suitable for circulation through closed systems.

Temperature range -10°C to 60°C

Stability: ±0.25°C or ±0.5°C (model dependent)

Choice of models with different cooling power - from 1300 to 3000W

Efficient, reliable and cost-effective alternative to cooling with mains water

Customised units with bespoke specifications also available



Lockable wheels allow RC units to be moved easily and ensure that they stay once put in position.

Digital controller for accurate and reproducible temperature setting. User-selectable high and low temperature alarms.

> Robust construction using corrosion resistant materials - long term durability and reliability in demanding applications.



Useful TUNE facility

enables automatic optimisation of the chillers closed-loop temperature control parameters to meet specific user requirements.

#### Inbuilt safety features

protects the user, equipment and application from over temperature, under temperature and flow failure.

### Applications

- Electronics cooling system for etch baths, glass coating for top-up display in aircrafts
- Industry print head cooling for textile industry, calibration system probe
- Academia physics and astronomy lab equipment cooling, sea water cooling for producing ikatite minerals
- Research seed research, cooling of scientific X-ray analytical units, SEM cooling

## Grant RC recirculating chillers

Technical specifications

		F	RC1400G	RC3000G
Dimensions	hxdxw mn	٦	655 x 93	36 x 483
Capacity		L	2.5	1.1
Temperature range	٥(		-10 t	0 60
Stability	±°C		0.25**	0.5***
Flow rate (max)	L/mir	۲	۱	5
Pump pressure (max)	mba	r	1.	6
Cooling power (typical)	@20°C V	/	1300	3000
	@0°C V	<b>/</b>	600	1500
	@-10°C V	<b>V</b>	150	575
Heater power	V	<i>√</i>	1500	_*
Overall consumption	220/240V V	<b>/</b>	3000	2000
Display			LE	ED
Display resolution	٥(		1.	.0
Electrical supply	N	/	230 (	50Hz)
Safety:				
temperature switchable under temperature the	rmostat			•
temperature fixed over temperature cut-out			•	-
level flow-fail device				
Refrigerant			RI	34a
EMC emissions	Clas	S	А	В
Weight	k	9	53	88

\* RC3000G has no heater so can only control against a heat load

\*\* With 10 litres of water in the system

\*\*\* With 25 litres of water in the system

#### Options and accessories

RC BYP	Bypass to overcome flow restrictions (flow <1 L/min), e.g. in narrow tubes or small cells.
RC PR	Pressure gauge to assist with setting up cooling systems and monitoring performance.
PRES	Priming reservoir to simplify priming in a closed loop system which has no filling port available on the RC inlet.
External probe	For remote sensing temperature control. On request only. Specify when ordering.
RC HF9	Rear connecting fittings (pair) 9mm internal diameter tube sizes respectively.
RC HF12	Rear connecting fittings (pair) 12mm internal diameter tube sizes respectively.
RC HF17	Rear connecting fittings (pair) 17mm internal diameter tube sizes respectively.

### Contact us today

Grant Instruments (Cambridge) Ltd w. www.grantinstruments.com 29 Station Road, Shepreth, Cambridgeshire, SG8 6GB

- t. +44 (0) 1763 260 811
- ► GrantInstruments
- y GrantInstrument
- e. salesdesk@grantinstruments.com in grant-instruments-cambridge-ltd