new art laboratories





SANLI Virus Specimen Collection Tube

The SANLI Medical transport medium is a complete solution specifically developed for the collection, transportation and storage of viral specimens – including SARS-CoV-2, influenza and enteroviruses.

- Can be stored at room temperature
- ✓ Effective antibiotic-antifungal combination enables various applications
- Conical transport tubes no decanting necessary for centrifugation
- 2 mL of transport medium
- Flocked swab for efficient oropharyngeal specimen collection

Rapid Tests


Laboratory Diagnostics


Laboratory Service

Consulting & Service





Test procedure











SANLI Medical Virus Transport Medium

Foetal bovine serum is added to a solution of Hanks' salts, maintaining virus activity over a wider temperature range. This facilitates the extraction of nucleic acid for subsequent samples, as well as the isolation of the virus.

Product code	Tube dimen- sions	Swab dimensions	Medium	Application	Product description	Pack size
780027	13x75 mm (7.5 mL)	Flocked swab Ø 5 mm, stick Ø 2.5 mm, length 150 mm	2 mL VTM	Oropharyn- geal	Conical tube packed with a flocked swab with predetermined breaking point for oropharyngeal specimen collection.	20 units / box

Our **customer service team** is on hand to answer any questions you may have by telephone, via email or in a personal consultation. For details on **your local sales office and regional free phone numbers, visit www.nal-vonminden.com.**



Rapid Tests

Laboratory Diagnostics

Laboratory Service

Consulting & Service

Rev02.01/XXC Rev02.01/074

Headquarters

 $\textbf{nal von minden GmbH} \cdot \textbf{Carl-Zeiss-Strasse} \ 12 \cdot 47445 \ Moers \cdot \textbf{Germany} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Roland Meißner} \cdot \textbf{Thomas Zander} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Roland Meißner} \cdot \textbf{Thomas Zander} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Roland Meißner} \cdot \textbf{Thomas Zander} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Roland Meißner} \cdot \textbf{Chief Executive Officers: Sandra von Minden} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Roland Meißner} \cdot \textbf{Chief Executive Officers: Sandra von Minden} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Roland Meißner} \cdot \textbf{Chief Executive Officers: Sandra von Minden} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Chief Executive Officers: Sandra von Minden} \\ \textbf{Chief Executive Officers: Sandra von Minden} \cdot \textbf{Chief Executive Officers: Sandra von Minden} \\ \textbf{Chief Executive Officers: Sandra von Minden Von Minden} \\ \textbf{Chief Executive Officers: Sandra von Minden Von Minden$

Commercial reg. Kleve · HRB 5679

Tel.: +49 941 290 10-0 · Fax: +49 941 290 10-50