

### new art laboratories



# nal von minden Alcohol tests

Breath, saliva and urine rapid tests Breathalysers for professional use





## nal von minden Alcohol tests

As the most ubiquitous and most underestimated drug there is, alcohol can also be considered the most dangerous. Alcohol can be an unexpected component in many everyday products. It is often unknowingly ingested during the consumption of certain foods, medications and even the use of cosmetics. Using the nal von minden alcohol tests, alcohol can be detected in breath, saliva and urine samples. Even measurements suitable for use by police can be obtained using our professional breathalyser, and our laboratory analyses deliver legally sound results.

In line with our motto 'We help you care!', we have the right test for you — whatever your needs may be.





### **Alcohol rapid tests**

#### Analysis of urine, saliva and breath samples

Road traffic control, custody battles, employment law, medical examinations, addiction therapy and transplantation medicine: Alcohol markers play a key role in various contexts, can even alter the course of a person's life. For this reason, a high degree of attentiveness is required in its detection. The nal von minden Drug-Screen<sup>®</sup> Alcohol rapid tests are your first step in this process.

#### **Advantages**

- Results within minutes
- ✓ Simple and hygienic test procedure
- ✓ Fast and reliable results
- Various multi-test combinations available for the simultaneous evaluation of numerous parameters
- In conformity with DIN EN ISO 15189 for medical laboratories, our rapid tests are available barcoded for use in reader-based evaluations
- Traffic control and customs

#### **Detection of alcohol (ALC)**

The direct detection of alcohol (ALC) is a reliable method of establishing alcohol consumption. Those who don't have a breathalyser to hand can make use of a fast and simple alcohol rapid test which is suitable for use with urine, saliva and breath (with or without balloon) samples. Our alcohol rapid tests detect acute intoxication — anything above 0.2 ‰. Depending on individual factors and the quantity consumed, alcohol can be detected in urine, saliva and breath samples for up to 12 hours.

#### Alcohol rapid test: Saliva or urine

| Description | Cut-off | Product code | Tests per pack |
|-------------|---------|--------------|----------------|
| Saliva test | 0.2 ‰   | 360006       | 50             |
| Urine test  | 0.2 ‰   | 101036       | 50             |

#### Alcohol screening breath test\*

| Description   | Product code | Tests per pack |
|---|--------------|----------------|
| Breath test with balloon 0.2 ‰  | 360010       | 25             |
| Breath test without balloon 0.2 ‰   | 360034       | 25             |
| Breath test without balloon 0.5 ‰   | 360035       | 25             |
| Breath test with balloon 0.5 ‰  | 360036       | 25             |
| *Coveral further versions and pack sizes are available. Den't besitete to contact vel |              |                |

\*Several further versions and pack sizes are available. Don't hesitate to contact us!

Have you heard about our EtG-ALC combi-tests for the simultaneous detection of ethyl glucuronide and alcohol, as well as the differentiation between current intoxication and recent consumption? Interested? Don't hesitate to get in touch!



#### **Advantages**

- Easy to use
- Fast results in just 5-10 minutes
- Minimum 11 month shelf life
- Can be integrated into multi-test combinations
- In conformity with DIN EN ISO 15189 for medical laboratories, our rapid tests are available barcoded for use in reader-based evaluations

#### **Drug-Screen EtG + ALC Test**

| Product                       | Format        | Sample | Cut-off                    | Product code | Pack size |
|-------------------------------|---------------|--------|----------------------------|--------------|-----------|
| Drug-Screen<br>EtG + ALC Test | Test cassette | Urine  | EtG 500 ng/mL,<br>ALC 0.2‰ | 102060       | 30        |



### **Detection of ethyl glucuronide (EtG)**

The nal von minden Drug-Screen<sup>®</sup> Ethyl glucuronide rapid test enables the qualitative detection of the direct, alcohol-specific biomarker EtG in human urine, with a cut-off of 500 ng/mL. In contrast to the detection of alcohol itself (which relies on its acute presence), or indirect biomarkers that only signal long-term excessive alcohol consumption, EtG covers the detection window in between. It is an excellent aid for the introduction or monitoring of treatment measures, but is also suitable for use in forensics. The EtG rapid test is an *in-vitro* diagnostic tool intended for professional use. Depending on environmental factors specific to the individual, EtG can be detected for an average of 1-4 days following the consumption of alcohol, which is longer than alcohol itself can be traced in urine. Our EtG rapid tests are available as test strips for immersion in urine samples, as well as test cassettes with pipettes provided. Both test formats are evaluated by comparing the colour intensity of the test line with the examples on the evaluation card.

#### **Advantages**

- Qualitative immunoassay
- Detection window: 1-4 days
- Easy to use
- Fast results in just 5-10 minutes
- Minimum 11 month shelf life
- In conformity with DIN EN ISO 15189 for medical laboratories, our rapid tests are available barcoded for use in reader-based evaluations

#### Note:

The evaluation cards are batch-specific and must be replaced with each new one.



#### EtG

Rapid test for the detection of ethyl glucuronide in urine

| Product         | Format        | Sample material | Cut-off   | Product code | Pack size |
|-----------------|---------------|-----------------|-----------|--------------|-----------|
| Drug Screen EtG | Test strips   | Urine           | 500 ng/mL | 101062       | 50        |
| Drug Screen EtG | Test cassette | Urine           | 500 ng/mL | 102062       | 30        |



### **Breathalysers**

Police officers in numerous countries throughout Europe use our breathalysers as part of traffic control. They are also suitable for use by justice authorities, clinics, pharmacies, companies, clubs and doctors' offices. Our alcohol breath tests are easy to use and of the highest quality, delivering reliable results in just a few seconds.





#### AlcoQuant 6020

This innovative breathalyser sets standards when it comes to safety. The electrochemical sensor targets alcohol, and short testing times mean that multiple tests can easily be carried out in succession. The use of the AlcoQuant 6020 by police authorities across the world is testament to its reliability.



#### **Advantages**

- Highly precise electrochemical sensor
- Simple, 3-button operation
- 🗸 Large, back-lit display
- Four-digit result display
- Streamlined design, lightweight and easy to handle
- Modern look
- Secure data transfer
- 🗸 Hygienic

#### Technical data:

- Measurement range: 0.00 to 5.50 ‰
- Accuracy: +/- 5 % at 1.00 ‰
- Evaluation time: Results within seconds
- Battery life: More than 1600 tests
- Dimensions: 190 mm x 61 mm x 38 mm
- Weight: 275 g incl. batteries
- Calibration: Recommended every 6-12 months
- No calibration constraint

#### Kit contents:

- 1x AlcoQuant 6020
- 25 mouthpieces
- 4 batteries
- 1 carry case
- User manual

| Product*            | Product code |
|---------------------|--------------|
| AlcoQuant 6020      | 360011       |
| 100 mouthpieces     | 360012       |
| Calibration service | 300200       |

\*Additional accessories and all pricing available upon request

#### **ALCOSCAN II (Basic plus)**

The ALCOSCAN II (Basic plus) is a digital breathalyser which uses the alcohol concentration in breath samples to calculate the blood alcohol concentration. With its high degree of accuracy, the electrochemical sensor is on par with professional devices and comes at an attractive price.



#### **Advantages**

- Highly precise electrochemical sensor
- Simple, 3-button operation
- Large, back-lit display
- Streamlined design, lightweight and easy to handle
- A cost-efficient alternative amongst devices with electrochemical sensors
- Maximum precision, short testing time
- 🗸 Storage card
- Warning signal for insufficient breath sample

#### Technical data:

- Measurement range: 0.00 to 4.00 ‰
- Accuracy: +/- 0.5 % at 1.0 ‰
- Evaluation time: 10 seconds at 1 ‰
- Battery life: 1500 tests
- Dimensions: 106 mm x 50 mm x 20 mm
- Weight: 100 g incl. battery
- Calibration: Every 6 months or after 500 tests

#### Kit contents:

- ALCOSCAN II breathalyser
- 5 mouthpieces
- 2 AAA batteries
- Belt bag
- Robust plastic carry case
- User manual in English & German

| Product*                            | Product code |
|-------------------------------------|--------------|
| ALCOSCAN II Basic plus breathalyser | 360014       |
| 50 mouthpieces                      | 360009       |
| Calibration service                 | 300200       |

\*Additional accessories and all pricing available upon request

#### Dräger Alcotest<sup>®</sup> 6820 med

The Dräger Alcotest<sup>®</sup> 6820 med is specially tailored for use in diagnostics. For example, it may be used for the diagnosis of vascular injuries following the application of an ethanol-based rinsing solution. It is also suitable for use in emergency situations or as part of methadone substitution programmes.



#### **Advantages**

- CE-certified in accordance with EU guideline 93/42/EWG for medical products
- Easy to use, short reaction time
- Minimal standard deviation
- Hygienic thanks to mouthpieces with a rebreathing lock, the subject's lips do not come into contact with the breathalyser
- Can be used even if subject is unconscious
- Increased storage capacity, logs up to 2000 test results

#### Technical data:

- Measurement range: 0 to 2.5 mg/L
- Accuracy: 0 to 0.5 mg/L: ± 0.008 mg/L
- Evaluation time: Results within a few seconds
- Dimensions: 65 mm x 147 mm x 39 mm
- Weight: 260 g incl. battery
- Calibration: Every 6 months or after 500 tests

#### Kit contents:

- ACE<sup>®</sup> ALP-1 med breathalyser
- 5 mouthpieces
- 2 AA batteries
- 1 CD-ROM with user manual MAX30 CD
- USB data cable
- Carry case
- Quickstart guide

| Product*                              | Product code |
|---------------------------------------|--------------|
| Dräger Alcotest <sup>®</sup> 6820 med | 300600       |

\*Additional accessories and all pricing available upon request

#### ACE ALP-1 med

The ACE ALP-1 med is also medically certified. As it is versatile and of high quality, it is suitable for use by both

authorities and those working in clinical settings such as hospitals, ambulances or doctors' practices.



#### **Advantages**

- CE-certified in accordance with EU guideline 93/42/EWG for medical products
- High quality, accuracy and versatility
- Comprehensive software menu, optional printer (available separately), PC software for test result evaluation
- Mouthpiece attachable at both ends

#### Technical data:

- Measurement range: 0.0-5.0 ‰
- Accuracy: ± 0.05 ‰
- Evaluation time: 10 seconds at 1 ‰
- Storage capacity: 500 results
- Dimensions: 133 x 64 x 33 mm
- Weight: 198 g incl. batteries
- Calibration: Every 6 months

#### Kit contents:

- ACE<sup>®</sup> ALP-1 med breathalyser
- 5 mouthpieces
- 2 AA batteries
- 1 CD-ROM with user manual MAX30 CD
- USB data cable
- Carry case
- Quickstart guide

#### Name\*

ACE® ALP-1 med

#### **Product code**

360064

\*Additional accessories and all pricing available upon request



Our range also contains further breathalysers. Get in touch with us on: +49 (0)941 29010-0





### **Laboratory Analysis**

Our nal von minden Drug-Screen<sup>®</sup> products are intended for use in pre-screenings. In situations with legal implications, or with potentially serious medical consequences, it is the norm for further laboratory analyses to be carried out.

#### How it works:

- 1. You can order an analysis through us online, or by telephone or fax.
- 2. We send you a sampling kit with all the information and materials you'll need for sending your sample to the lab.
- 3. Once the sample has been collected in accordance with the kit's instructions, you seal it in the pre-addressed envelope and send it to the laboratory.

In such cases, gas chromatography–mass spectrometry (GC-MS) or liquid chromatography–mass spectrometry (LC-MS) are the preferred methods. Should you require a laboratory analysis, we would be happy to arrange this for you.

- 4. Your results will be available within approximately 14 days, or within 5 working days for an express analysis.
- 5. The evaluation will be sent to you at the address provided.

| Description             | Product code |
|-------------------------|--------------|
| Laboratory analysis EtG | 640107       |



#### Germany

**Moers (Headquarters)** 

Carl-Zeiss-Strasse 12, 47445 Moers Tel.:....+49 2841 99 820-0 Fax: .....+49 2841 99 820-1

. . . . . . . . . . . . . . . .

#### Belgique / België

#### Français:



Fax gratuit: .....0800 747 07 dept.ventes@nal-vonminden.com www.nal-vonminden.fr Nederlands:

Fax: .....+31 70 3030 775 verkoop.dept@nal-vonminden.nl www.nal-vonminden.nl

Numéro gratuit:.....0800 718 82

#### Česká republika

czech@nal-vonminden.com www.nal-vonminden.cz

#### Danmark

Gratis nummer: ......80888753 Fax: .....+31 7030 30775 distribution@nal-vonminden.dk www.nal-vonminden.dk

#### Deutschland

Gratis-Fax: ......0800 589 2023 vertrieb@nal-vonminden.de www.nal-vonminden.de

#### España

Tel. directo: ......+49 941 290 10-759 ventas@nal-vonminden.es www.nal-vonminden.es

#### **Regensburg (Headquarters)**

| Friedenstrasse 32, 93053 Regensburg |
|-------------------------------------|
| Tel.:+49 941 290 10-0               |
| Fax:+49 941 290 10-50               |

#### France

Numéro vert: .....0800 915 240 Fax gratuit: .....0800 909 493 N° international: ..+49 941 290 10-739 Fax international: ...+49 941 290 10-50 dept.ventes@nal-vonminden.fr www.nal-vonminden.fr

#### Ireland

| www.nal-vonminden.com |              |
|-----------------------|--------------|
| sales-dept@nal-vonm   | inden.com    |
| Freefax:              | 1800 555 080 |
| Freephone:            | 1800 555 008 |

#### Italia

Numero verde: ...... 800 789 129 vendite.dipt@nal-vonminden.com www.nal-vonminden.it

#### Latinoamérica

Tel. directo:.....+49 941 290 10-761 ventas.latam@nvm-diagnostics.com www.nal-vonminden.es

#### Luxembourg

Fax gratuit: ......800 261 79 dept.ventes@nal-vonminden.com www.nal-vonminden.fr

#### Nederland

Gratis Tel: .....0800 022 2890 Gratis Fax: .....0800 024 95 19 verkoop.dept@nal-vonminden.nl www.nal-vonminden.nl

#### Norge

Free phone: ......800 16 731 norway@nal-vonminden.com www.nal-vonminden.com/no

#### Österreich

| www.nal-vonminden.de          |  |
|-------------------------------|--|
| vertrieb.abt@nal-vonminden.de |  |
| Gratis Fax:0800 298 197       |  |
| Gratis Tel:0800 291 565       |  |

#### Göttingen

| Robert-Bosch | -Breite 23, 37079 Göttingen |
|--------------|-----------------------------|
| Tel.:        | +49 551 270706-0            |
| Fax:         | +49 551 270706-50           |

#### Polska

| www.nal-vonminden.pl           |  |  |  |
|--------------------------------|--|--|--|
| sprzedaz.dept@nal-vonminden.pl |  |  |  |
| Bezpłatny fax:00800 4911594    |  |  |  |
| Bezpłatny tel.:00800 4911595   |  |  |  |

#### Portugal

| www.nal-vonminden.pt           |  |  |
|--------------------------------|--|--|
| vendas@nal-vonminden.pt        |  |  |
| Tel. direto:+49 941 290 10-735 |  |  |
| Tel. gratuito:                 |  |  |
|                                |  |  |

#### Schweiz / Suisse / Svizzera

#### Deutsch:

| Gratis-Tel.:                   |  |  |
|--------------------------------|--|--|
| Gratis-Fax:0800 837 476        |  |  |
| info@nal-vonminden.de          |  |  |
| www.nal-vonminden.de           |  |  |
| Français:                      |  |  |
| Numéro gratuit:0800 564 720    |  |  |
| Fax gratuit:0800 837 476       |  |  |
| dept.ventes@nal-vonminden.com  |  |  |
| www.nal-vonminden.fr           |  |  |
| Italiano:                      |  |  |
| Numero verde:0800 832 2146     |  |  |
| Fax gratuito:0800 837 476      |  |  |
| vendite.dipt@nal-vonminden.com |  |  |
| www.nal-vonminden.it           |  |  |

#### Suomi

| Ilmainen Puh:            | 0800 918263 |  |
|--------------------------|-------------|--|
| Ilmainen Fax:            | 0800 918262 |  |
| myynti@nal-vonminden.com |             |  |
| www.nal-vonminden.c      | om/fi       |  |

#### **Sverige**

sweden@nal-vonminden.com www.nal-vonminden.se

#### United Kingdom

Freephone: .....0808 234 1237 sales-dept@nal-vonminden.co.uk www.nal-vonminden.com

info@nal-vonminden.com www.nal-vonminden.com