

Integr^aline

HIGHLY RELIABLE ULT FREEZER WITH TRUE DUAL™ COOLING SYSTEM

Ultimate storage of your valuable samples with Arctiko TRUE DUAL™ cooling technology that will provide extra security for your samples.

Due to the independent cooling systems the freezer will maintain a temperature of at least -70°C, thereby samples will not be affected by unexpected failures.

- Independent **TRUE DUAL™** cooling system protecting the samples
- Reaching -80°C in less than 5 hours
- Direct cooling for optimized temperature stability and uniformity
- Low energy consumption made possible with ECO mode and ultra-efficient Vacuum Insulated Panels

- Adjustable shelves to accommodate any storage need you might have
- Easily removable inner doors make the defrost and maintenance procedure easy



MODEL	ULUF P390	ULUF P500	ULUF P610	ULUF P820
Item code	DAI 1800	DAI 1805	DAI 1810	DAI 1815
Capacity (L)	381	496	610	815
Exterior dimensions WxDxH (mm)	755x1002x1999*	895x1002x1999*	1035x1002x1999*	1285x1002x1999*
Interior dimensions WxDxH (mm)	466x620x1320	606x620x1320	746x620x1320	996x620x1320
Samples (pcs) in cryovials and 10x10 Cryoboxes	25.200	33.600	42.000	58.800

* Depth of 1002 mm includes detachable airflow management and noise reduction system

MODEL	ULUF P390	ULUF P500	ULUF P610	ULUF P820
DISCOUNT	25% Discount	25% Discount	25% Discount	25% Discount
Time frame for campaign	Until 30 June 2020	Until 30 June 2020	Until 30 June 2020	Until 30 June 2020

[1] Ready for immediate delivery

[2] Only applicable for orders accepted before June 30, 2020

[3] If the items are not in stock, we will assess the validity on the date of order and payment.

[4] We reserve the right to change specifications without notice. Subject to confirmation, availability and errors.

Contact for more information:

042-3009130, 08-400 200 42, 031-277001

E-Mail:

ORDER@LABTEAMET.COM