

04^{K30} Dry Bath Incubator



Specification

Model	K30
Temperature control range	RT.+5°C-100°C
Timer range	1min~99h59min
Temp. control accuracy	±0.5°C
Heating time (20°C to 100°C)	≤15min
Certificate	CE
Dimensions (L×W×H)mm	225×160×197
Weight	3kg

Blocks for K30



Ordering information

No.	Code	Description
01	AS-01030-00	K30, AC220V/AC120V, 50/60Hz, 100W

Accessories

No.	Code	Type	Description
01	AS-01031-01	A-30	96×0.2ml tubes
02	AS-01031-02	B-30	29×0.5ml+32×1.5ml tubes
03	AS-01031-03	C-30	61×0.5ml tubes
04	AS-01031-04	D-30	41×1.5ml tubes
05	AS-01031-05	E-30	41×2.0ml tubes
06	AS-01031-06	F-30	61×Φ6mm tubes
07	AS-01031-07	H-30	28×Φ12mm tubes
08	AS-01031-08	I-30	28×Φ13mm tubes
09	AS-01031-09	J-30	21×Φ15mm tubes
10	AS-01031-10	K-30	21×Φ16mm tubes
11	AS-01031-11	L-30	17×Φ18mm tubes
12	AS-01031-12	M-30	14×Φ20mm tubes
13	AS-01031-13	N-30	9×Φ25mm tubes
14	AS-01031-14	O-30	4×Φ32mm tubes