

06 MK200 Series High Temperature Dry Bath Incubator

Features

MK200 series high temperature dry bath incubator, the highest temperature of MK200-1 and MK200-2 can reach 150°C. The highest temperature of MK200-4 can reach 120°C



MK200-1



MK200-1A



Blocks for MK200 Series



MK200-2



MK200-4

Ordering information

No.	Code	Description
01	AS-01080-00	MK200-1, AC220V/AC120V, 50/60Hz, 200W
02	AS-01090-00	MK200-1A, AC220V/AC120V, 50/60Hz, 200W
03	AS-01100-00	MK200-2, AC220V/AC120V, 50/60Hz, 400W
04	AS-01110-00	MK200-4, AC220V/AC120V, 50/60Hz, 600W

Specification

Model	MK200-1	MK200-1A	MK200-2	MK200-4
Temperature control range	RT.+5°C-150°C	RT.+5°C-150°C	RT.+5°C-150°C	RT.+5°C-120°C
Timer range	1min-99h59min	1min-99h59min	1min-99h59min	1min-99h59min
Tem. stability @100°C	±0.5°C	±0.5°C	±0.5°C	±0.5°C
Temp. control accuracy @40°C	±0.3°C	±0.3°C	±0.3°C	±0.3°C
Temp. control accuracy @120°C	±1°C	±1°C	±1°C	±1°C
Heating time (40°C to 150°C)	≤30min	≤30min	≤30min	≤30min
Display accuracy	0.1°C	0.1°C	0.1°C	0.1°C
Dimensions (L×W×H)mm	280 × 240 × 110mm	280 × 240 × 110mm	280 × 240 × 110mm	380 × 240 × 110mm
Block quantity	1	1	2	4
Display screen	LED	LED	LED	LED

Accessories

No.	Code	Type	Description	Dimensions (mm)
01	AS-01081-01	BK01	49×Φ6mm	96×96×49mm
02	AS-01081-02	BK02	49×Φ7mm	96×96×49mm
03	AS-01081-03	BK03	25×Φ10mm	96×96×49mm
04	AS-01081-04	BK04	25×Φ12mm	96×96×49mm
05	AS-01081-05	BK05	25×Φ13mm	96×96×49mm
06	AS-01081-06	BK06	12×Φ15mm	96×96×49mm
07	AS-01081-07	BK07	16×Φ15mm	96×96×49mm
08	AS-01081-08	BK08	12×Φ16mm	96×96×49mm
09	AS-01081-09	BK09	16×Φ16mm	96×96×49mm
10	AS-01081-10	BK10	12×Φ19mm	96×96×49mm
11	AS-01081-11	BK11	16×Φ19mm	96×96×49mm
12	AS-01081-12	BK12	9×Φ20mm	96×96×49mm
13	AS-01081-13	BK13	9×Φ26mm	96×96×49mm
14	AS-01081-14	BK14	4×Φ28mm	96×96×49mm
15	AS-01081-15	BK15	4×Φ40mm	96×96×49mm
16	AS-01081-16	BK16	49×0.5ml tubes	96×96×49mm
17	AS-01081-17	BK17	25×1.5ml tubes	96×96×49mm
18	AS-01081-18	BK18	25×2.0ml tubes	96×96×49mm
19	AS-01091-01	BK19	96×0.2ml tubes	78×114×26mm
20	AS-01091-02	BK20	96×0.2ml ELISA plate	81×123×19mm

BK19 and BK20 can only be used in MK200-1A