

04 MX100-4A Microplate Shaker

Description

The MX100-4A is designed for shaking ELISA plate(96 well plates), TCPS and 1.5/2.0ml centrifuge tubes

Features

1. Continuous or timing operation
2. DC brushless motor for gentle and powerful vibration
3. Orbital of 3mm for efficient mixing in small sample vessels
4. Simple spring holders for easy and quick plate installation



MX100-4A

Specification

Model	MX100-4A
Speed Setting Resolution	10rpm
Speed	200~1500rpm @Standard / 200~1000rpm @Deep well plate or tube holder
Timer range	1min~99h59min
Mixing Orbit	3mm
Input Power	AC100V -240V
Dimensions (mm)	290 x 270 x 120 mm
Weight	6.5kg

Ordering information

No.	Code	Description
01	AS-04040-00	MX100-4A Microplate shaker
02	AS-02041-01	20x1.5ml tube holder (apolegamy)
03	AS-02041-02	20x0.5ml tube holder (apolegamy)