

# 02 Mini-15K High Speed Centrifuge

## Features

1. Extremely quiet operation and easy-to-follow LCD display
2. Capacity: 12 x 1.5/2.0ml tubes
3. Low sample heating (only 12°C after 30min, at max.speed)
4. Autoclavable rotor (121°C, 20min), and includes standard rotor made of anodized aluminum
5. RPM/RCF setting as required
6. Timer can be set up to 99min
7. The small size of the unit makes it convenient for benchtop use or in a hood or in refrigerated environments



**01**  
Operation interface



**02**  
Rotor



Mini-15K

## Specification

Model	Mini-15K
Max. rotational speed	14,500rpm
Max. RCF	14,000 x g
Timer	20sec-99min
Acceleration time to max. speed	15s
Braking time from max. speed	15s
Dimensions (D×W×H)mm	270 × 230 × 155mm
Power	100W

## Ordering information

Code	Description
AS-08110-00	Mini-15K High speed centrifuge, 50/60Hz, 3A

## Accessories

AS-08041-01	Adapters for 0.5 ml tubes
AS-08041-02	Adapters for 0.2 ml tubes