

01 MiniT Series Biological Indicator Incubator



MiniT-1

Description

The MiniT-1 and MiniT-3 Biological Indicator Incubator are designed for temperature incubation of the EO Biological Indicators. Both of them are compact, tabletop incubator.

Features

1. LCD display
2. User can set the incubating temperature and time.
3. Exact correspondence of the set and actual temperature
4. It is easy to set up and use
5. The Incubator has many wells for placement of the biological indicator vials and the wells can be custom-made.



MiniT-3

Specification

Model	MiniT-1	MiniT-3
Temp. Control Range	RT.+5°C~80°C	RT.+5°C~100°C
Timer Range	1min-99h59min	1min-99h59min
Temp. Control Accuracy	±0.5°C	±0.5°C
Display Accuracy	±0.1°C	±0.1°C
Heating Time(25 to 57°C)	≤6min	≤6min
Dimensions (mm)	110 x 150 x 80	120 x 152 x 112
Weight	0.5kg	0.85kg

Ordering information

No.	Code	Description	
01	AS-09020-00	MiniT-1, AC100~240V, 50/60Hz, 0.5A	
02	AS-09030-00	MiniT-3, AC100~240V, 50/60Hz, 0.5A	
Accessories			
No.	Code	Type	Description
01	AS-09021-01	A-MiniT-1	15 x Φ10mm
02	AS-09021-02	B-MiniT-1	15 x Φ12mm
03	AS-09021-03	C-MiniT-1	6 x Φ15mm
04	AS-09021-04	D-MiniT-1	6 x Φ16mm
05	AS-09021-05	E-MiniT-1	24 x Φ8.8mm
06	AS-09021-06	F-MiniT-1	18 x Φ9mm
07	AS-09021-07	G-MiniT-1	12 x Φ12.5mm
08	AS-09021-08	H-MiniT-1	6 x Φ16.5mm