


01 OS-100/OS-200 Orbital Shaker

Features

The OS-100 and OS-200 Orbital Shaker are powerful variable speed shaker which provides efficient orbital motion. Speed and time are under microprocessor control. The OS-200/OS-100 shaker can be applicable in different labs: microbiology, chemistry, immunology, biochemistry and molecular biology



PP-1




PP06-250



PP12-100



UP-1



MP-1

01 Platforms for OS-100 / OS-200

Specification

| | OS-100 | OS-200 |
|-------------------|--------------------|--------------------|
| Model | OS-100 | OS-200 |
| Shaking Speed | 50--300rpm | 50--250rpm |
| Orbit | 10mm | 20mm |
| Timer range | 1min~99h59min | 1min~99h59min |
| Load capacity | See below table | See below table |
| Max. capacity | 3.0kg | 2.5kg |
| Motor | DC brushless | DC brushless |
| Operating ambient | 4-45°C | 4-45°C |
| Dimensions (mm) | 280 x 270 x 110 mm | 280 x 270 x 110 mm |

Ordering information

| No. | Code | Description |
|-----|-------------|--|
| 01 | AS-04030-00 | OS-100, AC100-240V, 50/60HZ, 2A |
| 02 | AS-04020-00 | OS-200, AC100-240V, 50/60HZ, 2A |
| 03 | AS-04021-01 | PP-1, flat platform with non-slip rubber mat |
| 04 | AS-04021-02 | PP12-100, platform with clamps for flasks, 100-150ml |
| 05 | AS-04021-03 | PP06-250, platform with clamps for flasks, 250-300ml |
| 06 | AS-04021-04 | UP-1, universal platform with adjustable bars for different types of flasks, bottles |
| 07 | AS-04021-05 | MP-1, platform with springs |