

Equipment For Gel Card

K37-24 Gel-Card Incubator

■ Features

1. Incubator for 24-Gel Cards in racks at 37°C
2. 2 independent incubation zones with a capacity for 12 Gel cards
3. Digital interface screen for control and visualization of remaining incubation time and real temperature for each zone
4. Fixed pre-set temperature of 37 degree
5. Pre-set incubation time of 15 minutes, which can be adjusted by the user
6. Audible programmable alarm to notify end of incubation period

■ Specification

Capacity	24 Gel cards
Temperature	37°C
Timer	15min (Variable)
Power	AC220V or AC120V 50~60Hz, 200W
Dimensions	380 x 280 x 140 mm (D x W x H)

Code	Description
AS-10020-00	Gel-Card Incubator AC220V or AC120V 50~60Hz 3A



K37-24

TD2-12/TD2-24 Gel-Card Centrifuge

■ Features

1. Silent centrifuge for up to 24 or 12 ID-cards.
2. Removable and exchangeable spinning head
3. Automatic balance control
4. Pre-defined centrifugal parameters for optimum Gel-Cards processing; acceleration, speed, braking and time
5. Digital interface visor control and visualization of centrifuge speed and time remaining for processing
6. Programmable audible alarms for end of centrifuge time periods

■ Specification

Model	TD2-12	TD2-24
Capacity	12 Gel Cards	24 Gel Cards
Speed Range	800~1800rpm	500~1500 rpm
Centrifugal Force	54~272g	33~294g
Timer	1~20min	1~20min
Power	AC100~240V, 85W	AC100~240V, 85W
Dimensions (DxWxH)	440 x 340 x 170mm	530 x 410 x 180mm
Net Weight	7kg	11.5kg

Code	Description
AS-10010-00	Gel-Card Centrifuge AC100V~240V 50~60Hz 3A
AS-10030-00	Gel-Card Centrifuge AC100V~240V 50~60Hz 3A



TD2-12



TD2-24



Gel-Cards

