

MX IIR



Ear Loop Medical Face Mask Type II R

The Haika medical face mask offers excellent protection to limit the transmission of infection.

The 3 ply mask consists of three layers of non-woven fabric including an anti-bacterial outer layer, one meltblown middle layer and a soft inner layer.



QTY: 50



NON-STERILE | TYPE IIR

Key Features and Benefits

Excellent Bacterial Filtration Efficiency 98%

Splash Resistant

Easy to fit and adjusted over the nose and mouth

Low breathing resistance and excellent comfort to wearer

Tested in accordance of EN14683: 2019

Ordering Information:

- **Item No:** HAK8001
- **Packaging:** 50 per Dispenser Box
40 Dispensers Per Case

Product Specification

- **Size:** 17.5cm x 9.5cm
- **Material:** 3ply - Ear loop

Compliant With:

- Medical Device Directive
- MDD 93/42/ECC
- 2007/47/EC
- EN 14683: 2019
- CLASS 1

www.haikagloves.com



W W W . l a b t e a m e t . c o m

Labutrustning • Service • Kvalitet

LabTeam Scandinavia AB
Adress Vasatorpsvägen 1
SE-254 57 Helsingborg Sverige

Telefon
Direkt 042-300 91 32
Kundtjänst 042-300 91 30
Stockholm 08-400 200 42
Göteborg 031-277 001
Danmark 36 98 02 30

Mail labteamet@labteamet.com
Order order@labteamet.com
Service service@labteamet.com

Web
www.labteamet.com