

NX476



Nitrile Examination Gloves

Non-latex, powder free, non-sterile

HAIKA™ NX476 are lightweight, ambidextrous, single-use medical nitrile examination gloves.

They offer healthcare professionals high performance barrier protection as well as supreme comfort and touch sensitivity.

Key Features

Secure Control

A textured surface enhances slip resistance and allows for a secure grip.

Supreme Comfort and Sensitivity

The lightweight and elastic nature of the nitrile delivers an excellent sense of touch and a comfortable wearing experience even when worn for long periods.

Effective Barrier Protection

Helps to protect users and patients from the transmission of infections.

Beaded Cuff

Durable beaded cuff design enables the glove to be donned and removed efficiently and securely.

Powder-free and Latex-free

Eliminates powder-induced irritation and makes it ideal for people with a Type I latex sensitivity.

Minimises the Risk of Allergic Contact Dermatitis

Free of Thiurams, Mercaptobenzotiazoles, Thiazoles and Thioureas.



EU Compliance

- MDR (EU) 2017/745 Class I
- PPE Regulation (EU) 2016/425 CAT I
- EN455

Recommended For

- Accident and Emergency
- Bio-science
- Dental
- Laboratory
- Pharmacy
- Wards

NX476



Nitrile Examination Gloves

TECHNICAL

- **Type:** Non-sterile, powder-free single use examination gloves
- **Material:** Nitrile Butadiene Rubber (NBR)
- **Colour:** Light Blue
- **Glove Shape:** Flat with a beaded cuff
- **Impermeability:** AQL 1.5
(AS PER EN 455-1)
- **Force at Break (N):** >6 N
(AS PER EN 455-2)
- **Shelf Life:** 5 years
(IN ORIGINAL PACKAGING WHEN STORED AS PER DIN 7716, ISO 2230)
- **Country of Origin:** Malaysia

Material Information

- **Premium Performance Nitrile Butadiene Rubber**
- **Powder-free**
- **Beaded cuff**
- **Light Blue**

Sizing

- **Overall Length Size:** 240mm
(AS PER EN455-2)
- **Wall Thickness:** Min. 0.09 mm
(MEASURED IN SINGLE LAYER OF THE PALM AREA)

Ordering Information

- **Product Code:** HAK4761 S
HAK4762 M
HAK4763 L
- **Pack Quantity:** 100 single gloves by weight, in dispensers
- **Case Quantity:** 10 dispensers per case

General Information

Conformance:

HAIKA™ NX468 medical examination gloves have been tested and found to comply with all parts of EN455

Registering Authority: In Europe the gloves are CE marked indicating compliance with the PPE regulation 2016/425 (CAT 1)

They are also in compliance with Medical Directive 93/42/EEC Medical Device Class 1 (non-sterile)



Labutrustning • Service • Kvalitet

LabTeam Scandinavia AB
Adress Vasatorpsvägen 1
SE-254 57 Helsingborg Sverige

Telefon
Direkt 042-300 91 32
Kundtjänst 042-300 91 30
Stockholm 08-400 200 42
Göteborg 031-277 001
Danmark 36 98 02 30

Mail labteamet@labteamet.com
Order order@labteamet.com
Service service@labteamet.com

Web
www.labteamet.com