

CSPP Controller

MOTORIZED PIPETTE CONTROLLER

Features:

Light and fatigue free: weighs only 167g including the battery

• Lithium battery with no memory effect and up to 12 hours of continuous use

• Can be charged while in use

Powerful: 50mL pipette is filled in less than 10 seconds

• Express dispensing: depress the dispensing button

beyond the built-in point of resistance

Gravity delivery and blow-out dispensing mode selection

Fully adjustable speed-setting

• Exchangeable 0.2 μm PTFE disc filter acts as effective shut-down system if overfill occurs

• Exchangeable pipette holders are made of silicon

• For use with glass or plastic pipettes ranging from 0.1mL to 100mL

• Vapour outlet for venting out corrosive liquids or fumes

Fully autoclavable nose cone

UV-resistant housing

Complies with CE and RoHS



Ordering Information

Cat. No. Description

PA-100 CAPPController Motorized Controller, 0.1-100mL