

CAPP Solo

MULTICHANNEL PIPETTES

CAPP Solo multichannel pipettes offer industry-standard pipetting accuracy and precision in a very affordable package. Manufactured with budget-constrained customers in mind, these pipettes offer quick, reproducible pipetting without compromising on quality.

CAPP Solo multi-channel pipettes are supplied in 8-channel and 12-channel variants that are guaranteed to last long in academic, diagnostic and other comparable laboratories that require workhorses that can withstand extended usage cycles.

Features:

- Colour-coded for quick identification of multi-channel pipettes with different volume capacities
- Remarkably lightweight pipette chassis for its class for more comfortable pipette handling
- The precise counter mechanism offers accurate volume setting and can be locked in place to avoid accidental turning while pipetting
- Accurate liquid dispensing across all channels for repeatable results
 - Durable build quality ensures chemical resistance, bump protection, and thermal stability
- Long-term calibration retention even after sterilising in an autoclave
 - Each tip cone is universally compatible with tips from the most common brands
 - Supplied with a standard two-year warranty



Ordering Information

CAPP Solo Multichannel Pipettes

Cat. No.	Description
S10-8	CAPP Solo 8-channel pipette, 0.5 – 10 μ L
S100-8	CAPP Solo 8-channel pipette, 10 – 100 μ L
S300-8	CAPP Solo 8-channel pipette, 30 – 300 μ L
S10-12	CAPP Solo 12-channel pipette, 0.5 – 10 μ L
S100-12	CAPP Solo 12-channel pipette, 10 – 100 μ L
S300-12	CAPP Solo 12-channel pipette, 30 – 100 μ L

