

The **CAPP** *Tempo* pipette controller has revolutionary features that separate it from other pipette controllers.

LED Display

• Display shows aspiration and dispense speeds

- Displays remaining battery life and flashes when the power is low
- Display turns-off after 20 seconds of non-use to conserve battery

Adjustable speeds

- Powerful motor fills 25mL in 5 seconds (at speed level six)
- Separate speed controls for aspiration and dispense modes
- Six different speed settings for both aspiration and dispense modes

Electronic Brake System (EBS)

- Apply EBS to slow down at the last second to dial in the most accurate volume
- Press EBS a second time to return to max speed
- EBS maximizes speed and precision

Power Supply

- Lithium ion battery, full charge provides 8 hours of continuous use
- Battery can be charged while pipette is in use

Ordering Information

| Cat. No. | Description |
|----------|--------------------------|
| Т100-В | CAPPTempo Blue |
| T100-G | CAPPTempo Green |
| T100-R | CAPPTempo Red |
| T100-V | CAPPTempo Violet |
| T100-CS | CAPPTempo Charging Stand |

www.capp.dk - +45 6613 6140 - info@capp.dk

the the the