

# EcoSpin ES-2 / ES-4



## Key features

The EcoSpin™ Centrifuge has especially been developed to centrifuge EcoPlates™. The instrument is easy to handle and due to its clever design, doesn't need any maintenance or calibration. Once delivered, it is ready to use.

Speed: 3,000 rpm (constant)  
Interval: 12 sec. (constant)  
Environmental Temperature: 4°C to 50°C  
Relative Humidity at 20°C: 80 %  
Power: 24 V, 1 A, max. 24 Watt  
Size: 90 x 170 x 120 mm  
Weight: 2,5 kg

**Optional rotors to  
operate 2 or 4  
plates at a time!**

## Interface

- Two optional rotors (2 plates / 4 plates)
- Mechanical lid lock
- Super fast start/stop
- Preprogrammed optimal speed and time of a spindown
- Automatically starts when closing the lid
- Emergency lid unblocking mechanism
- Main power switch on/off



For more information please see a "how to" guide:  
[www.youtube.com/elmimedia/videos](http://www.youtube.com/elmimedia/videos)