

# Fit of Eppendorf Storage Boxes to Eppendorf ULT Freezers and Freezer Racks

Eppendorf Storage Box				Eppendorf Freezer Racks											
Height	Grid	Order No.	Intended vessels	2 inch racks				3 inch racks				4 inch racks			
				Article number / boxes per rack				Article number / boxes per rack				Article number / boxes per rack			
				K0641-1900 / 16	K0641-1690 / 13	K0641-3003 / 10	K0641-3000 / 28	K0641-1890 / 12	K0641-1700 / 9	K0641-3004 / 6	K0641-3001 / 20	K0641-1880 / 8	K0641-1750 / 6	K0641-3005 / 4	K0641-3002 / 12
2 inch	10 x 10	0030 140.508	cryogenic tubes	■	■	■	■	■	■	■	■	■	■	■	■
	9 x 9	0030 140.516	cryogenic tubes	■	■	■	■	■	■	■	■	■	■	■	■
	8 x 8	0030 140.524	1-2 mL tubes	■	■	■	■	■	■	■	■	■	■	■	■
2.5 inch	5 x 5	0030 140.532	5 mL (snap cap) tubes	-	-	-	-	■	■	■	■	■	■	■	■
3 inch	9 x 9	0030 140.540	cryogenic tubes	-	-	-	-	-	■	■	■	■	■	■	■
	5 x 5	0030 140.613	5 mL (screw cap) tubes	-	-	-	-	-	■	■	■	■	■	■	■
4 inch	9 x 9	0030 140.567	cryogenic tubes	-	-	-	-	-	-	-	-	■	■	■	■
5 inch	3 x 3	0030 140.583	15 mL conical tubes	-	-	-	-	-	-	-	-	-	-	■	■
		0030 140.591	50 mL + (4 x 15 mL) conical tubes	-	-	-	-	-	-	-	-	-	-	■	■

The outer footprint dimensions are identical for all boxes: 128,9 x 128,9 mm

## Eppendorf ULT Freezer Fitting Legend

CryoCube® F570/h, Premium U410, HEF® 410

CryoCube FC660/h, Innova® C585 and C760

Innova U101

Innova U360, U535 and U725/-G

■: Yes

