## **Pressure tank 80**

Cat. no: ZC-0080



## Destiny:

Metal pressure tank used for water storage demineralized water.

## **Building:**

The tank is built of a fully corrosion resistant housing made of fiberglass. The construction of the tank guarantees no contact of stored water with external environment (light, air).

- Height: 78 cm
- Diameter: 39 cm
- Total volume: 88,5 l
- Total capacity at pressure 2,8 bar: 76,2 l
- Max. pres. full collection: 0,35 MPa
- Ambient temperature: 2 40°C
- Empty tank pressure: 0,4-0,5 bar

