

# Control Plus

## Upright Control Plus Solid Door Refrigerator 353L

GDPA code **PPSR353**



 **MADE IN BRITAIN**

### Key features

- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms



#### New controller

New updated Controller system for improved user experience. Read, Record, Reset simply with the touch of one button, with no stress of accidentally resetting on the clear digital display.



#### Dual probe

Greater temperature accuracy given from two temperature probe readings. One probe measuring air temperature and one measuring load (or product) temperature.



#### Porthole

Factory fitted 16mm porthole enables quick and easy access for any secondary temperature monitoring probes.



#### Flexible internal storage

Multiple shelf locations along with a wire basket provides a flexible storage solution.



#### Anti-microbial handle

Powder coated to help prevent cross contamination by inhibiting the growth of microbes.

## Specification: Upright Control Plus Solid Door Refrigerator 353L

Gross Capacity (litres)	353L
Net Capacity (litres)	349L
Unit External Width (mm)	595mm
Unit External Depth (mm)	640mm
Unit External Height (mm)	1860mm
Product Weight (kg)	69kg
Packed Width (mm)	650mm
Packed Depth (mm)	727mm
Packed Height (mm)	1885mm
Packed Weight (kg)	75kg
Country of Origin	UK
Internal Dimensions (mm) (WxHxD)	480mm W x 1682mm H x 445mm D
Temperature Range	2°C to 8°C (tolerance + or - 2)
Colour	White/green
External Finish	Painted steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted steel
Door Type	Solid
Door Style	Hinged, reversible
Door Liner	Flat
Door Open Depth (mm)	1185mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet (number & position)	2 adjustable front feet, 2 rear rollers
Porthole	Yes
Porthole Size	16 mm
Porthole Location	Rear Bottom
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	6
Number of Shelf Positions	14
Shelf Size (mm) (WxDxH)	385mm W x 6mm H x 500mm D
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Basket Size (mm) (WxDxH)	
Interior Light	Yes
Light Type	N/A
Light Activation	N/A
Fan	Yes
Spark Free	No
Defrost setting	Auto
Controller Type	Digital, external
Number of Probes	2
Probe Type	Air & Load (Glycol Solution)
Ambient Temperature Range	10 to 32
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	220-240
Insulation Thickness (mm)	60mm
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Wall Brackets	No
Optional Drawer System	No
Optional Digital Lock	No
Other Door Option	Glass PPGR353UK - 444410665