





# MHCare Medical Product Catalog

Please contact us for special prices  
[info@myhealthcare.com.tr](mailto:info@myhealthcare.com.tr)

**LAMINATED COVERALL (EU) 2017/745 - 93/42EEC (non-sterile)****Ref No: 470801****Contents:**


- Zipper front
- Available sizes: S~3XL
- Elastic wrists, ankles and waist for a good seal and comfortable fit
- Low linting fabric
- Protect skin and clothing against aerosols, light spray, splashes, chemicals and dust
- Laminated material for resistance to liquid, abrasion and tearing
- Uses: Pharmaceutical industries, food industries, inspections, coverage of patients in diagnostic examinations, guided tour.

**COVERALL (EU) 2017/745 - 93/42EEC (non-sterile)****Ref No: 470301****Contents:**

- Zipper front
- Available sizes: S~3XL
- Elastic wrists, ankles and waist for a good seal and comfortable fit
- Low linting fabric
- Protect skin and clothing against aerosols, light spray, splashes, chemicals and dust
- Uses: Pharmaceutical industries, food industries, inspections, coverage of patients in diagnostic examinations, guided tour.

**BOUFFANT (non-sterile)****Ref No: MHCAP01-B1****Contents:**


- Suitable to contain even thick hair
- Equipped with a vertical rear opening
- Easy adaptability to the garment
- Breathable
- Material: Spunbond + Spunbond (Non-woven) material
- One size

**HEADCAP (WITH CIRCULAR ELASTIC)  
(non-sterile)****Ref No: MHCAP01-CR****Contents:**


- Suitable to contain even thick hair
- Equipped with a vertical rear opening
- Easy adaptability to the garment
- Breathable
- It ensures adequate and constant sealing without causing irritation
- The seams resistant to traction and tearing
- Anti-reflective colors
- Minimum weight : 18 GSM
- One size

**DRUM SHAPE HEADCAP (non-sterile)****Ref No: MHCAP01-DR****Contents:**


- Suitable to contain even thick hair
- Elastic behind the neck
- Breathable
- It ensures adequate and constant sealing
- Anti-reflective colors
- The seams are resistant to traction and tearing
- Minimum weight : 25 GSM
- One size

**HEADCAP (FOR MAN) (non-sterile)****Ref No: MHCAP01-MAN****Contents:**


- Suitable to contain even thick hair
- Laces behind the neck (sufficient length to ensure easy and suitable fastening)
- Equipped with a vertical rear opening
- Easy adaptability to the garment
- Breathable
- It ensures adequate and constant sealing
- Anti-reflective colors
- Minimum weight : 25 GSM
- One size

**INTEGRAL HEADCAP (non-sterile)****Ref No: MHCAPO1-IT****Contents:**


- Suitable to contain even thick hair
- Full-length model headgear
- Perforated print to ensure maximum breathability
- Cover the head, forehead, ears up to the neck with wide bands
- Equipped with a comfortable back fastening
- Equipped with a lacing that completely encloses the hair
- Anti-reflective colors
- Minimum weight : 25 GSM
- One size

**LAMINATED INTEGRAL HEADCAP (non-sterile)****MHCAPO1-LIT****Contents:**

- Suitable to contain even thick hair
- Full-length model headgear
- Suitable to cover head and shoulders
- Protect skin and clothing against aerosols, light spray, splashes, chemicals and dust
- Laminated material for resistance to liquid, abrasion and tearing
- Material: SS+ PE (Spunbond + Spunbond +Polyethylene)
- Dimension: 40 x 48cm ( Hole inside: 22cm diameter)
- Minimum weight : 56 GSM
- One size

**SHOE COVER (non-sterile)****Ref No: MHCOVER01-CSC-PE****Contents:**


- Provided with non-slip sole
- Ample fit
- Waterproof
- Overshoes with elastic at the ankle
- High tear resistance
- Material:PE
- Minimum foot length : 45 CM
- One size

**SHOE COVER (ANKLE-HEIGHT)  
(non-sterile)****Ref No: MHCOVER01-SC-M****Contents:**

- Provided with non-slip sole
- Ample fit
- Waterproof
- In TNT and PE / PP
- Overshoes with elastic at the ankle
- High tear resistance
- Minimum foot length : 45 CM
- One size

**SHOE COVER CLASSIC (non-sterile)****Ref No: MHCOVER01-CSC****Contents:**


- Provided with non-slip sole
- Ample fit
- Overshoes with elastic at the ankle
- Material: Spunbond + Spunbond (Non-woven) material
- Minimum foot length : 45 CM
- One size

**SHOE COVER (KNEE-HEIGHT)  
(non-sterile)****Ref No: MHCOVER01-SC-L****Contents:**

- Non-woven fabric
- Breathable
- Knee-length with a wide fit
- It has an elastic closure and reinforced central seam
- High tear resistance after prolonged use
- Minimum foot length : 45 CM
- One size

**SCRUB SUIT (SHORT SLEEVE)**  
(non-sterile)

Ref No: MHGOWN01 SCR

**Contents:**

- Medical blue color for SMS fabric
- Pattern Type: Solid Color
- Collar: V-neck
- Sleeve length: Short Sleeve
- Gender: Unisex
- Size: S~3XL
- Weight: 35GSM
- Material: SMS

**FILTERING HALF MASK FFP2 NR (With 10 different color options)**

Ref No: MHMSK-01

**Contents:**

With 10 different color options

- Adjustable nose clip
- Ultrasonic welded straps
- Soft elastic rope
- Implementation standard:  
EN149:2001+A1:2009 FFP2 NR /  
CE1463

**DISPOSABLE PATIENT GOWN**  
(non-sterile)

Ref No: MHGOWN01-PT


**Contents:**

- Water repellent (for sizes M, L, XL)
- Equipped with rear lacing, with neck and trocar sleeve closure system
- Non-Woven fabric
- Minimum weight : 30 GSM
- Size: XS~3XL



### DISPOSABLE STANDARD GOWN (non-sterile)

Ref No: MHGOWN01




#### Contents:

- Water repellent (for sizes M, L, XL)
- Equipped with back lacing and neck closure system
- Equipped with long sleeves with elasticized cotton cuff
- Non-Woven fabric
- Minimum weight : 30 GSM
- Size: XS~3XL

### DISPOSABLE SANITARY STAFF GOWN (sterile)

Ref No: MHGOWN01-ST




#### Contents:

- Water repellent (for sizes M, L, XL)
- Equipped with back lacing and neck closure system
- Equipped with long sleeves with elasticized cotton cuff
- Non-Woven fabric
- Minimum weight : 30 GSM
- Size: XS~3XL


### DISPOSABLE SANITARY STAFF GOWN (W/Reinforced) (non-sterile)

Ref No: MHGOWN01-ST-RI




#### Contents:


- PE / PP Double coupled gown (laminated) (for sizes M, L, XL)
- Totally waterproof
- Equipped with back lacing and with a neck closure system
- It ensures adequate barrier against blood and liquids when subjected to pressure
- Non-Woven fabric
- Minimum weight : 40 GSM
- Size: XS~3XL

**DISPOSABLE SANITARY STAFF GOWN  
(W/Reinforced) (sterile)****Ref No: MHGOWN01-ST-RIS****Contents:**

- PE / PP Double coupled gown (laminated) (for sizes M, L, XL)
- Waterproof with back fastening and collar closure system
- Equipped with long sleeves with stretch cotton cuff
- Non-Woven fabric
- Minimum weight : 40 GSM
- Size: XS~3XL

**VISITOR GOWN (non-sterile)****Ref No: MHGOWN01-VS****Contents:**



- Water repellent (for sizes M, L, XL)
- Equipped with back lacing and with collar closure system
- Long sleeves with elastic or cuff
- Non-Woven fabric
- Minimum weight : 17 GSM
- Size: XS~3XL

**STANDARD SURGICAL GOWN  
(non-sterile)****Ref No: MHGOWN01 SG****Contents:**

- SMS fabric construction
- Equipped with back lacing and neck closure system
- Equipped with long sleeves with elasticized cotton cuff
- Non-Woven fabric
- Minimum weight: 30 GSM
- Size: XS~3XL

**Full Reinforced SURGICAL GOWN  
(non-sterile)**

Ref No: MHGOWN01 FRI





**Contents:**

- Medical blue color for SMS fabric
- Critical areas reinforced with Polyethylene
- Elastic cuffs
- Gender: Unisex
- Disposable
- Latex free
- Size: S-XXL
- Weight: 35-40GSM

**Reinforced SURGICAL GOWN  
(non-sterile)**

Ref No: MHGOWN01 RI



**Contents:**

- Medical blue color for SMS fabric
- Critical areas reinforced with Polyethylene
- Elastic cuffs
- Gender: Unisex
- Disposable
- Latex free
- Size: S-XXL
- Weight: 35-40GSM



## ABOUT US

MHCare Medical LTD. is a fully licensed Medical Device supply company proudly based in Turkey. We pride ourselves in being one of the largest medical personal protective equipment suppliers in Turkey, Europe and North America.

We specialize in high-quality medical personal protective equipment which aligns with our mission to protect our end-users from harmful substances and infectious diseases!

Nothing is more important to us than providing a perfect service for our customers. We have consistent levels of inventory on hand ensuring immediate turnaround times for our clients' PPE orders

We strive to ensure on-time delivery of orders 100% of the time with a strong and diverse network of logistics partners, we can handle all logistics requirements by air, sea, rail, or ground

Our employees have business, biopharma/biochemistry, QA/QC degrees, bachelor of arts, science degrees, and more. We have the capacity to customize products to suit specific end user requirements

 /myhealthcareturkey

 /myhealthcareturkey

 /company/myhealthcareturkey

Değirmiçem District 16023 Street No: 5A Şehitkamil/Gaziantep, Turkey

[www.myhealthcare.com.tr](http://www.myhealthcare.com.tr) - [info@myhealthcare.com.tr](mailto:info@myhealthcare.com.tr)

 +90 501 337 61 27