





FILTERING HALF MASK
FFP2 NR



FFP2 mask with
different color options



EN149:2001+A1:2009



FFP2 NR


MELTBLOWN
PROTECTIVE MASK

With 10 different color options

Features



MHMSK-01
EN149:2001+A1:2009
FFP2 NR CE 1463



- Adjustable nose clip
- Ultrasonic welded straps
- Soft elastic rope
- 5 layer filtration, double layer melt-blown spray filtration. (DOP \geq 95% PFE \geq 95% BFE \geq 99%)



FFP2 NR

MELTBLOWN
PROTECTIVE MASK

Implementation standard:
EN149:2001+A1:2009 FFP2 NR

CE Certificate NO:
CW/PPER/8/04/2021

Monolithic packaging of FFP2 NR

1pc/bag:

- 5 languages on the package (English, German, Italian, Spanish and French)
- Transparent package design
- QR code available for checking certification and product detail



Option 1	10pcs/box	680pcs/ctn
Option 2	25pcs/box	900pcs/ctn

Carton size: 60x40x40cm

Truck	
With 24 Pallets	Option 1/408000pcs
	Option 2/540000pcs
Without Pallet	Option 1/540000pcs
	Option 2/720000pcs

* Each pallet includes 25ctn



Monolithic packaging of FFP2 NR

10pc/bag:

- 5 languages on the package (English, German, Italian, Spanish and French)
- QR code available for checking certification and product detail



10pcs/box - 680pcs/ctn

Truck

With 24 Pallets	408000pcs
Without Pallet	540000pcs

Carton size: 60x40x40cm

* Each pallet includes 25ctn



Monolithic packaging of FFP2 NR

25pc/bag:

- 5 languages on the package (English, German, Italian, Spanish and French)
- QR code available for checking certification and product detail



25pcs/box - 900pcs/ctn

Carton size: 60x40x40cm

Truck


With 24 Pallets	540000pcs
Without Pallet	720000pcs

* Each pallet includes 25ctn

Certificates



CE 1463 Certificate (Module B)



CE 1463 Certificate (Module C2)



ISO 13485 Certificate



Scan for
more details



ABOUT US

MHCare Medical LTD. is a fully licensed Medical Device supply company proudly based in Turkey. We pride ourselves in being one of the largest medical personal protective equipment suppliers in Turkey, Europe and North America.

We specialize in high-quality medical personal protective equipment which aligns with our mission to protect our end-users from harmful substances and infectious diseases!

Nothing is more important to us than providing a perfect service for our customers. We have consistent levels of inventory on hand ensuring immediate turnaround times for our clients' PPE orders

We strive to ensure on-time delivery of orders 100% of the time with a strong and diverse network of logistics partners, we can handle all logistics requirements by air, sea, rail, or ground

Our employees have business, biopharma/biochemistry, QA/QC degrees, bachelor of arts, science degrees, and more. We have the capacity to customize products to suit specific end user requirements



/myhealthcareturkey



/myhealthcareturkey



/company/myhealthcareturkey

Değirmiçem District 16023 Street No: 5A Şehitkamil/Gaziantep, Turkey

www.myhealthcare.com.tr - info@myhealthcare.com.tr

+90 501 337 61 27