# new art laboratories





### **Flocked Swab**

For the collection of nasal and nasopharyngeal specimens. The flocked swab can be used for several test methods – including PCR, rapid tests and many more.

### Your advantages at a glance:

- Specially designed for the collection of nasal and nasopharyngeal specimens
- Multifunctional use for several test methods
   (e.g. SARS-CoV-2 PCR, influenza/RSV rapid tests, etc.)
- Flocked with nylon
- ✓ Predetermined breaking point at 85 mm
- Individually packed and sterile (CE0197-certified)

Rapid Tests

Laboratory Diagnostics

Laboratory Service

Consulting & Service

## new art laboratories




#### **Flocked Swab**

#### Product details:

Product code	Swab type	Swab dimensions	Application	Pack size
780026	Flocked nylon	Tip Ø 2.80 mm, stick Ø 0.80 mm - 2.35 mm, length 150 mm, breaking point at 85 mm	nasal and nasopharyngeal	1 pc

**Original size** 



#### Do you know about our VTM solution?

With the CITOSWAB® Virus Collection and Transport Kit, nal von minden also offers a complete solution - including swab, transport medium and transport tube - for transportation and preservation of viral specimens, such as COVID-19, influenza, swine flu, RSV or adeno-viral based infections.

Our customer service team is on hand to answer any questions you may have by telephone, via email or in a personal consultation. For details on your local sales office and regional free phone numbers, visit www.nal-vonminden.com.



Commercial reg. Kleve  $\cdot$  HRB 5679 Tel.: +49 941 290 10-0 · Fax: +49 941 290 10-50