

iStir Uno / iStir Quattro

Motorless Ultra Slim Magnetic Stirrer



The motorless magnetic stirrer is a truly maintenance free ultra flat stirring device for various laboratory applications.

- ◆ No moving parts - truly maintenance free
- ◆ Low aspect ratio (just 12 mm thick)
- ◆ Software driven gradual acceleration
- ◆ Adjustable speed 15 to 1500 RPM (Uno) & 15 to 1200 RPM (Quattro)
- ◆ Microprocessor controlled with last run memory feature
- ◆ Pulse mode - reverses the direction every 30 secs for proper mixing
- ◆ Universal power supply - small wall mount adaptor
- ◆ IP 64 complaint design (Water spill resistant)
- ◆ Corrosion resistant SS top surface to fight harsh chemicals
- ◆ Stirring cap. 800 ml for Uno & 400 ml per station for Quattro

| Specifications | iStir Uno | iStir Quattro |
|-----------------------------|---|-------------------------------------|
| Max RPM | 1500 rpm | 1200 rpm |
| Speed Setting | 15 to 1500 RPM (at every 50 RPM) | 15 to 1200 RPM (at every 50 RPM) |
| Speed indicator | Yes | |
| Load capacity | 800 ml | 400 ml |
| Recommended Stir bar length | 25 mm | |
| Pulse function | Yes | |
| Increment | 50 RPM \pm on each press | |
| Last run memory | Yes | |
| Dimensions (W X D X H) | 120 x 164 x 12 mm | 457 x 126 x 13.4 mm |
| Weight | 400 gm | 1.9 kg |
| Protection Class | IP 65 compliant design (Water spill resistant) | |




**Extremely
slim
design**



**Pulse mode
for thorough
mixing**



**SS top for
corrosion
resistance**



**Microprocessor
controlled**

*OEM request can be entertained. Specifications can be changed without notice for quality improvement.