

## IDC/F CE



### DESCRIPTION

Infectious disease control shield (spit shield). This is a complete face shield assembly with elastic attachment and foam forehead cushion. Comes with 5 shields in a resealable bag.

---

### SPECIFICATION

<b>MATERIAL</b>	Mylar (Anti-Static)
<b>WEIGHT</b>	28 g (1 oz)
<b>COLOR</b>	Clear
<b>STANDARD</b>	OSHA 1910.1030, EN 166

---

### ORDER DATA

<b>PRODUCT NAME</b>	IDC/F CE
<b>PRODUCT NUMBER</b>	6256411

---

### APPLICATION

Medical, Dental, Corrections, Police, Fire, EMT, Industrial, Health & Beauty, Service & Hospitality, Education.

---

### COMPETITIVE ADVANTAGE

Protects against splash, particulate and bloodborne pathogens.

Anti-fog coating.

Each unit can be disposed of after use at minimal expense.

Re-sealable bag keeps the unused product clean and ready for use.

Individual units may be rolled into a tube and stowed in a cargo pocket for easy access when needed.

One size fits all.

# Product Data Sheet



Declaration of Conformity / Instruction for Use

DoC

IFU