

Data sheet

Cooled Incubator (IL) ILW 115 Smart PRO



The photo above is for reference only, may show additional options not included in standard equipment. The real appearance, particularly color and structure of the material may differ from the ones presented in the photo.

Advantages of the SMART PRO controller: - large (7"), clear, full colour touch screen

- LAN, USB ports and WiFi for communication and data transfer
- multi-segment time and temperature programs
- overview of data in tabular and graphic form
- visual and sound alarm
- Admin function for management
- password protected log-in
- internal memory for programs and data storage
- event registry with user notifications
- LabDesk software and user manual for direct download



Smart PRO - preview screen



Bank: ING BŚ o/Wodzisław Śl. Nr rachunku: 90 1050 1403 1000 0022 0879 5498 Nr rej. W BDO 000005609


NIP: 647-216-26-87 REGON: 276111984 Sąd Rejonowy w Gliwicach, X Wydział Gospodarczy KRS 0000083696



TECHNICAL DATA

air convection	forced	
chamber capacity [I]	112	
working capacity [I]	112	
controller	microprocessor PID	
display	7" full colour touch screen	
TEMPERATURE		
temperature range [°C]	0+100 / -10+100 (option)	
temperature resolution every [°C]	0,1	
temperature fluctuation at 37°C [+/-°C]*	0,2	
temperature variation at 37°C [+/-°C]*	0,3	
over temperature protection	class 3.3 to DIN 12880	
CHAMBER		
door type	double(5) / door with viewing window (option)	
interior		
Smart PRO	DIN 1.4301	
IG Smart PRO	DIN 1.4301	
housing		
Smart PRO	powder coated sheet	
IG Smart PRO	stainless steel linen finish	

vidth A	660
eight B	1140
epth C	720
nternal dims [mm]	
idth D	460
eight E	540
epth F	450



shelves (standard max)	2 7
max shelf workload [kg] /2/	25
- reinforced shelf version (PW) [kg] /3/	50
max unit workload [kg]	60
- reinforced unit version (W) [kg] /4/	120
weight [kg]	90





ELECTRICAL PARAMETERS

voltage**	230V 50-60Hz
nominal power [W]	500
refrigerant	1234ze / GWP=1
warranty	24 months
manufacturer	POL-EKO-APARATURA

all the above technical data refer to standard units (without optional accessories)

* - fluctuation measured in centre of the chamber; in space, variation (K) calculated for chamber as:

K= +/- (T average max. - T average min.) / 2

** - other power supplies on request

1 - depth doesn't include 50 mm of power cable

2 - on uniformly loaded surface

3 - reinforced shelf

4 - additional internal glass door

OPTIONS AND ACCESSORIES

Order number: */A INOX Order number: */P INOX Order number: */PP Order number: */PW Order number: KUW GN*/* Order number: QLK* Order number: */S or */S INOX Order number: */ST or */S INOX 潫 Order number: *PLUS

...

6

Order number: */T

door with viewing window

Stainless steel wire shelf INOX

Perforated shelf

Reinforced shelf

Stainless steel cuvettes

Castors

Table with castors

Base on castors

Automatic defrosting function

Low temperature version



Bank: ING BŚ o/Wodzisław Śl. Nr rachunku: 90 1050 1403 1000 0022 0879 5498 Nr rej. W BDO 000005609

NIP: 647-216-26-87 REGON: 276111984 Sąd Rejonowy w Gliwicach, X Wydział Gospodarczy KRS 0000083696



POL-EKO-APARATURA SP. J. A. POLOK – KOWALSKA, S. KOWALSKI 44-300 WODZISŁAW ŚLĄSKI, UL. KOKOSZYCKA 172 C



-

Order number: */W

Order number: OWW/OWW LED



+



Order number: PT 100



Order number: GNZ

Order number: */FIT

>0<



ALARM

USB

LAN

Order number: BPP 12

Order number: BRT/*/L or IQ/OQ/PQ

Order number: PORT ALARM

Order number: USBK

Order number: LANK

Order number: OCZ/20

Order number: OCZ/30

Order number: OCZ/60

20mm







(((●)))

Order number: OCZ/100

Order number: KD





Order number: CCTV CAM

Order number: CCTV

Reinforced version

Interior lighting

Phytotron system

Additional temperature sensor

Internal socket

Calibration and IQ, OQ, PQ qualification

Battery backup for display

Dry alarm contact

USB cable

LAN cable

Non-standard access port 20 mm

Additional access port 30 mm

Non-standard access port 60 mm

Non-standard access port 100 mm

Access control

Camera inside the chamber

Additional camera

NIP: 647-216-26-87 REGON: 276111984 Sąd Rejonowy w Gliwicach, X Wydział Gospodarczy KRS 0000083696



Bank: ING BŚ o/Wodzisław Śl. Nr rachunku: 90 1050 1403 1000 0022 0879 5498 Nr rej. W BDO 000005609