

Data sheet

Cooled Incubator (ST) ST 6 Smart



The photo above is for reference only, may show additional options not included in standard equipment. The real appearance, particularly color and structure of the material may differ from the ones presented in the photo.

Advantages of the SMART controller:

- 4,3", clear, full colour touch screen
- LAN, USB ports for data transfer
- multi-segment time and temperature programs
- visual and sound alarm
- internal memory for programs and data storage
- event registry
- user manual for direct download
- Quick change of program parameters
- Alarm Bar
- operating with gloves on



Smart - preview screen



Bank: ING BŚ o/Wodzisław Śl. Nr rachunku: 90 1050 1403 1000 0022 0879 5498 Nr rej. W BDO 000005609

NIP: 647-216-26-87 REGON: 276111984 Sąd Rejonowy w Gliwicach, X Wydział Gospodarczy KRS 0000083696




TECHNICAL DATA

air convection	forced
chamber capacity [l]	400
working capacity [I]	324
controller	microprocessor PID
display	4,3" full colour touch screen
TEMPERATURE	
temperature range [°C]	+3+40 / +70 (opcja)
temperature resolution every [°C]	0,1
temperature fluctuation at 37°C [+/-°C]*	0,3
temperature variation at 37°C [+/-°C]*	0,6
over temperature protection	class 1.0 to DIN 12880 / class 3.3 (option)
CHAMBER	
door type	solid / glass or double (option) /4/
interior	
B Smart	aluminium
C Smart	stainless steel to DIN 1.4016
CS Smart	stainless steel to DIN 1.4016
P Smart	DIN 1.4301
PS Smart	DIN 1.4301
housing	
B Smart	powder coated sheet
C Smart	powder coated sheet
CS Smart	stainless steel polished
P Smart	powder coated sheet
PS Smart	stainless steel polished

overall dims [mm] /1/

width A	620
height B	1900
depth C	650
internal dims [mm]	
width D	480
width D'	520
height E	1660
depth F	420
depth F'	480
depth G	320
height H	1440



shelves (standard max)	4 10
max shelf workload [kg] /2/	10
max unit workload [kg]	60
weight [kg]	90



Bank: ING BŚ o/Wodzisław Śl. Nr rachunku: 90 1050 1403 1000 0022 0879 5498 Nr rej. W BDO 000005609 NIP: 647-216-26-87 REGON: 276111984 Sąd Rejonowy w Gliwicach, X Wydział Gospodarczy KRS 0000083696



ELECTRICAL PARAMETERS

voltage**	230V 50-60Hz
nominal power [W]	350
refrigerant	R1234ze / GWP=1
warranty	24 months
manufacturer	POL-EKO-APARATURA

all the above technical data refer to standard units (without optional accessories)

* - fluctuation measured in centre of the chamber; in space, variation (K) calculated for chamber as:

K= +/- (T average max. - T average min.) / 2

** - other power supplies on request

1 - depth doesn't include 50 mm of power cable

2 - on uniformly loaded surface

3 - reinforced shelf

4 - additional internal glass door

OPTIONS AND ACCESSORIES

	Order number: */C	Internal glass door
	Order number: */A	External glass door
	Order number: */P	Wire shelf
	Order number: */P INOX	Stainless steel wire shelf INOX
	Order number: */PP	Perforated shelf
	Order number: */PW	Reinforced shelf
\bigcirc	Order number: KUW GN*/*	Stainless steel cuvettes
ALU	Order number: ST/CHL/SWP ALU	Aluminum drawer with powder coated slides
	Order number: ST/CHL/SWP INOX	Stainless steel drawer with powder coated slides
	Order number: ST/CHL/SWPN INOX	Stainless steel drawer with stainless steel slides



Bank: ING BŚ o/Wodzisław Śl. Nr rachunku: 90 1050 1403 1000 0022 0879 5498 Nr rej. W BDO 000005609 NIP: 647-216-26-87 REGON: 276111984 Sąd Rejonowy w Gliwicach, X Wydział Gospodarczy KRS 0000083696



	Order number: ORG-FARM	Pharma organizer
Ø	Order number: QLK*	Castors
ו••	Order number: *PLUS	Automatic defrosting function
+70°C	Order number: ST/70	Extended temperature range to 70°C
	Order number: */FOT	Photoperiodic system
	Order number: GNZ	Internal socket
LAB DESK	Order numberLabDesk	LabDesk software
>0<	Order number: BRT/*/L or IQ/OQ/PQ	Calibration and IQ, OQ, PQ qualification
DIN 3.3	Order number: */3.3	Over temperature protection 3.3 class according to DIN 12880
•	Order number: BPP 12	Battery backup for display
	Order number: PORT ALARM	Dry alarm contact
USB	Order number: USBK	USB cable
LAN	Order number: LANK	LAN cable
20mm	Order number: OCZ/20	Non-standard access port 20 mm
30mm	Order number: OCZ/30	Additional access port 30 mm
60mm	Order number: OCZ/60	Non-standard access port 60 mm
()) 100mm	Order number: OCZ/100 Bank: ING BŚ o/Wod:	
	Nr rachunku: 90 1050 Nr rej. W BDO 00000	1403 1000 0022 0879 5498 Sąd Rejonowy w Cliwicach, X Wydział Gospodarczy KRS 0000083696





Order number: KD

Access control

