

TECHNICAL DATASHEET

Deep well plates 1,2 ml - OV-7696 547

Features:	
Description	96 Deep well plate
Volume	1,2 ml
Dimension	See technical drawing Bottom shape: round (U)
Material	Polypropylen (PP) 100% virgin polymer, clean room class 8 production
Sterilization	Non-sterile
Purpose	Sample storage
Temperature range	- 80 °C to 121 °C (autoclavable)
Purity	Free from DNase/Rnase Free from human DNA Free from PCR inhibitors Free from endotoxins Free from ATP
Storage	Dry, at room temperature
Shelf life	3 years after month of production

Packaging


Sub packaging	4 pcs
Pieces per pack	32
Pack per EPAL	60

29.05.2020

TECHNICAL DRAWING

Deep well plates 1,2 ml - OV-7696 547

Deep well plates 1,2 ml - OV-7696 547



Dimensions : Millimeters

Revision 1: 29.05.2020