

# 23% MORE SAMPLES WITH ARCTIKO CRYOVIALS AND BOXES

Arctiko Cryovials and boxes are designed for storing biological material, human or animal cells at temperatures as low as  $-150^{\circ}\text{C}$

## ARCTIKO CRYOVIALS


- Slim design with external thread invented for storage in 10x10 boxes
- Certified Rnase-free, Dnase-free, Pyrogen-free, and DNA-free
- Sterilized after ISO 11137-2:2013 "Sterility Assurance Level"(SAL) by standard products  $10^{-6}$
- Star-foot grip for easy removal and fitting of the cryovial lid
- Uniquely designed lids - only  $\frac{1}{2}$  turn to fasten
- The additional storage capacity creates a direct reduction in CO2 discharge
- Arctiko Cryovials are available in 5 different sizes: 0.5ml, 1.0ml, 2.0ml, 4.0ml, 5.0ml



**ARCTIKO**  
Cryovials

## ARCTIKO CRYO BOXES 10X10

- The 10x10 boxes increase the storage capacity with 23% compared to conventional boxes
- Click lock in the lid ensures secure handling and samples will be stored optimally at all times
- Cryo boxes with numerical codes for quicker identification
- Standard measurements 130 x 130mm.
- Arctiko Cryo Boxes are available in 5 different heights: 32mm, 40mm, 50mm, 80mm, 90 mm



**ARCTIKO**  
Cryo Boxes

**ARCTIKO**  
true cooling specialists