




FREESTANDING FRIDGE-FREEZER COMBI SOLID DOOR 324L WHITE



- ✓ **Operating Temperature Range of +2°C to +10°C for Fridge and -18°C to -25°C for Freezer**
Conforms to MHRA guidelines
- ✓ **NEW Advanced Controller**
Records temperatures accurate to 0.1°C
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **2 Fitted Locks as Standard (one on fridge, one on freezer)**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**
For ease of movement
- ✓ **Wire Shelves**
Enables optimum air flow inside unit
- ✓ **Plastic Baskets and Door Organisers**
Maximises space and offers protection to contents

SKU CODE: PRODUCT CODE:

444441957 LSC324EU

DESCRIPTION:

Freestanding Fridge-freezer Combi Solid Door 324L White



FULL SPECIFICATION – LSC324EU FREESTANDING LABORATORY FRIDGE-FREEZER COMBI SOLID DOOR 324L WHITE

| | |
|---|--|
| Description | Freestanding |
| SKU Code | 444441957 |
| EAN Code | 5052263019577 |
| GDPA Code | LSC324EU |
| Gross Capacity (L) | 324 (Fridge: 190, Freezer: 134) |
| Net Capacity (L) | 302 |
| Unit External Width (mm) | 600 |
| Unit External Depth (mm) | 600 |
| Unit External Height (mm) | 1870 |
| Product Weight (kg) | 70 |
| Packed Width (mm) | 657 |
| Packed Depth (mm) | 734 |
| Packed Height (mm) | 1885 |
| Packed Weight (kg) | 78 |
| Country of Origin | UK |
| Internal Dimensions (mm) (WxDxH) | Fridge: 510 x 430 x 805, Freezer: 463 x 445 x 780 (step:250d x 260h) |
| External Finish | White |
| Internal Finish | White |
| Type of Plastic | SRL800 (HIPS/ABS Mix) |
| Temperature Range | Fridge: +2°C to +10°C, Freezer: -18°C to -25 °C |
| Colour | White/Green |
| Door Type | Solid |
| Door Style | Hinged, Reversible |
| Door Open Depth (mm) | 1190 |
| Replacement Door Seal | No |
| Open Door Sensor | Yes |
| Lock | Yes |
| Adjustable Feet | 2 adjustable feet, 2 rear rollers |
| Temperature Display on Panel | Yes |
| Temperature Alarm | Yes |
| Number of Shelves | 3 (Fridge) |
| Number of Shelf Positions | 6 |
| Shelf Size | 3: 320 x 505 x 4.5 |
| Shelf Material | Wire |
| Moulded Drawers | 2 Solid half (Fridge) |
| Adjustable Door Bins | 2 full + 4 half |
| Interior Light | No |
| Fan | No |
| Spark Free | Yes |
| Auto or Manual defrost | Fridge: Automatic, Freezer: Manual |
| Controller type | Digital, External |
| Climate Class | SN |
| Refrigerant Type | R600a |
| CFC Free | Yes |
| HFC Free | Yes |
| HCFC Free | Yes |
| EU Plug | Yes |
| Other Plug Types Available | Yes |
| Voltage Supply | 230 |
| Aux Output | Yes |
| Insulation Thickness (mm) | 60 |
| Insulation Type | Cyclopentane |
| Anti Microbial Handle | Yes |