

# Pharmacy Plus Bluetooth Enabled Refrigerator

## 400L Pharmacy Plus Upright Fridge with Glass Door with Left Hand Hinge


GDMA code      PPGR400BT-LHH

SKU code        444411318



### Key features

- ✓ Local monitoring and alarm history via Bluetooth
- ✓ Power failure alarm with backup battery powered recording
- ✓ Solid and glass door options
- ✓ Intelligent fan management system for enhanced performance



#### Bluetooth Connected Controller

Local temperature monitoring and alarm history via Bluetooth connectivity. Power failure alarm with battery back-up, providing continuous temperature monitoring in the event of power outage.

#### Lec Medical Smartphone App

Contactless local monitoring of air and load temperatures via smart phone app. Downloadable temperature report and alarm history.

#### Intelligent Fan Management

Smart fan optimisation aids in temperature stability following door openings, minimising fan activity when door is open to reduce the rate of temperature increase.

#### Dual Probe Temperature Monitoring

Our dual probe system offers greater temperature accuracy within your appliance, simultaneously measuring both air and product temperatures.

#### Flexible Internal Storage

Multiple shelf locations provides flexible storage solutions.

Specification: 400L Pharmacy Plus Upright Fridge with Glass Door with Left Hand Hinge

|                                     |  |
|-------------------------------------|--|
| Gross Capacity (litres)             | 394L                                   |
| Net Capacity (litres)               | 389L                                   |
| Unit External Width (mm)            | 595mm                                  |
| Unit External Depth (mm)            | 660mm                                  |
| Unit External Height (mm)           | 1935mm                                 |
| Product Weight (kg)                 | 69.6Kg                                 |
| Packed Width (mm)                   | 640mm                                  |
| Packed Depth (mm)                   | 660mm                                  |
| Packed Height (mm)                  | 1995mm                                 |
| Packed Weight (kg)                  | 78.6Kg                                 |
| Country of Origin                   | UK                                     |
| Internal Dimensions (mm) (WxDxH)    | 515mm x 480mm x 1783mm                 |
| Temperature Range                   | +2°C to +8°C                           |
| Colour                              | White/green                            |
| External Finish                     | Painted Steel                          |
| Internal Finish                     | HIPS*                                  |
| Type of plastic                     | SRL800 (HIPS/ABS mix)                  |
| Door Finish                         | PVC / Glass                            |
| Door Type                           | Glass                                  |
| Door Style                          | Hinged, Left hand only, Non-reversible |
| Door Liner                          | Flat                                   |
| Number of Door Trays                | n/a                                    |
| Number of Door Tray Positions       | n/a                                    |
| Door Open Depth (mm)                | 1150mm                                 |
| Replacement Door Seal               | Yes                                    |
| Open Door Sensor                    | Yes                                    |
| Lock                                | Yes, 2 Keys                            |
| Adjustable Feet (number & position) | 2 adjustable (front)                   |
| Porthole                            | Yes                                    |
| Porthole Size                       | 16 mm                                  |
| Porthole Location                   | Rear Bottom                            |
| Temperature Display on Panel        | Yes                                    |
| Temperature Alarm                   | Yes                                    |
| Remote Alarm Contacts               | Yes                                    |
| Audible Alarms                      | Yes                                    |
| Visual Alarms                       | Yes                                    |
| Number of Shelves/Drawers           | 7                                      |
| Number of Shelf Positions           | 12                                     |
| Shelf Size (mm) (WxDxH)             | 496mm x 405mm x 6mm                    |
| Shelf Material                      | Polythene Coated Wire                  |
| Drawer Size (mm) (WxDxH)            | n/a                                    |
| Drawer Material                     | n/a                                    |
| Number of Baskets                   | 1                                      |
| Basket Material                     | Polythene Coated Wire                  |
| Basket Size (mm) (WxDxH)            | 485mm x 255mm x 171mm                  |
| Interior Light                      | Yes                                    |
| Light Type                          | LED                                    |
| Light Activation                    | Door Switch / Controller               |
| Fan                                 | Yes - Tangential                       |
| Spark Free                          | no                                     |
| Defrost setting                     | Auto (Off Cycle)                       |
| Controller Type                     | Digital BT, External                   |
| Battery Back Up/Hours               | no                                     |
| Bluetooth                           | Yes                                    |
| NFC                                 | Yes                                    |
| WIFI                                | no                                     |
| Number of Probes                    | 2                                      |
| Probe Type                          | Air & Load (Silicone Oil)              |
| Ambient Temperature Range           | 10°C to 32°C                           |
| Refrigerant Type                    | R600a                                  |
| CFC Free                            | Yes                                    |
| HFC Free                            | Yes                                    |
| HCFC Free                           | Yes                                    |
| Ammonia Free                        | Yes                                    |
| UK 3 Pin Plug                       | Yes                                    |
| Other Plug Types Available          | Yes                                    |
| Voltage Supply                      | 220-240V                               |
| Insulation Thickness (mm)           | 45mm                                   |
| Insulation Type                     | High Density Polyurethane Foam         |
| Blowing Agent                       | Cyclopentane                           |
| Anti Microbial Handle               | Yes                                    |
| Wall Brackets                       | No                                     |
| Optional Drawer System              | No                                     |
| Optional Digital Lock               | No                                     |
| Other Door Option                   | No                                     |